INSTRUCTIONS FOR USE

REUSABLE PRESSURE CABLES

NOTE: Not all symbols must be shown on the label

CE Mark

Consult Instructions for Use or Electronic Instructions for Use

Non-Sterile

Manufacturer

Country of Manufacture
Date of Manufacture: Year-Month-Day

Catalogue Number

Lot Number

Keep Dry

Keep Away from Sunlight

Temperature Limit

Medical Device

Unique Device Identifier

Latex Free

Does not Contain Phthalates

Quantity

Authorized Representative in Switzerland

Authorized Representative in the European Union

Importer

CAUTION: US Federal law restricts this device to sale by or on the order of a physician

Separate collection as waste electrical and electronic equipment (WEEE)

Compatibility

Monitor Device Family

Sensor Device Family


Product Length

REUSABLE PRESSURE CABLES

This document applies to all Reusable Pressure Cables of Nissha Medical Technologies marketed under the brand name "Integral Process".

PRODUCT DESCRIPTION

1. Sensor-Side Connector
2. Cable
3. Monitor-Side Connector

Depending on the specific monitoring device or sensor, different cable types are available. The cables differ in monitor-side connector and sensor-side connector, cable material, and cable length. Pressure Cables are reusable, non-sterile, and non-invasive medical devices. They are suitable for multiple reuse with multiple patients when correctly processed according to these instructions for use.

INTENDED USE

Reusable Pressure Cables are accessories for invasive blood pressure (IBP) monitoring. The cables are used to connect an invasive blood pressure sensor (transducer) to an appropriate monitoring device.

INDICATIONS

Any medical condition requiring invasive patient blood pressure monitoring and/or diagnosis.

CONTRAINDICATIONS

There are no known contraindications for use.

INTENDED USER GROUP

Products may be used by medical professionals only who must have the necessary qualifications and training.

INTENDED PATIENT POPULATION

There are no patient population restrictions associated with the products.

WARNINGS AND PRECAUTIONARY MEASURES

- Carefully read these instructions for use before using the device. Not following these instructions could lead to harm to patients and/or users.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- Follow the instructions for use of the monitor/sensor being used.
- Any modification of the product by the end user renders the product liability of Nissha Medical Technologies null and void.
- Before use, ensure that the cable is used with a compatible monitor/sensor.
- Before use, visually check that the cable is in working order. There must be no broken or frayed wires or damaged connectors, and no cracks or discoloration. If damage is found, replace the cable.
- The product is not suitable for use in an MRI environment. Doing so may cause patient or user harm.

USE

- When connecting, disconnecting or handling the cable, always hold it by the connectors rather than the cable to avoid damage.
- Connect the cable to the sensor using the respective connector.
- Carefully route the cable to:
 - avoid the risk of patient entanglement,
 - avoid tripping hazards,
 - avoid direct contact between the cable and exposed skin.
- Connect the cable to the monitor using the respective connector.
- After use, clean the cable according to the processing instructions.

DISPOSAL

The Cables fall into the category of Waste from Electrical and Electronic Equipment (WEEE) and must be disposed accordingly. Please note that disposal must be in line with the applicable laws and requirements of the respective state.

CONDITIONS FOR OPERATION, TRANSPORT, AND STORAGE

	OPERATION	TRANSPORT	STORAGE
Temperature	10–40 °C / 50–104 °F	-40–70 °C / -40–158 °F	5–50 °C / 41–122 °F
Relative Humidity	30 % to 85 % noncondensing		

- Store the product in cool and dry place.
- Keep the product away from direct sunlight and from the immediate vicinity of light sources, heating, or other heat sources.

REGULATORY INFORMATION

- Report any serious incidents involving this product to the manufacturer and the competent authority in which the user and/or patient is established.

COMPATIBILITY

Further compatibility information can be found on the product website: <https://www.hs.nisshamedical.com/IFU>

PROCESSING

These instructions have been validated by the manufacturer as being capable of preparing the product for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

WARNINGS	• Do not allow any liquid or moisture to get on any of the connection points of the connectors. • Use only recommended cleaning and disinfecting substances referenced in this document. Others may reduce product lifetime, damage the cables or cause safety hazards.
LIMITATIONS ON PROCESSING	• The number of processing cycles shall be limited to no more than 300. • Cables should be discarded if any cable strain reliefs (sections between connectors and cable material) are soiled.
PREPARATION BEFORE CLEANING	• Visually check that the cable is in working order. If damage is found, dispose the cable.

CLEANING	• Clean the cable by wiping with a lint-free cloth moistened with a mild dishwashing detergent solution. For the solution, mix 1 teaspoon (5 milliliters) mild dishwashing detergent in 1 US gallon (3.8 liters) of water. • Do not allow any liquid or moisture to get on any of the connection points of the connectors. Visually verify that any soiling, stains or residue are removed. • Repeat cleaning if soiling, stains or residue are still observed. • Remove cleaning residue by wiping with a lint-free cloth moistened with distilled water. • Wipe dry with a clean lint-free cloth.
DISINFECTION	• Follow the disinfectant manufacturer's preparation instructions as applicable. • Fill a dip tank with a 2% to 3% hydrogen peroxide disinfectant (such as Revital-Ox® RESERT® High Level Disinfectant or ANIOXYDE 1000). • Partially immerse the cable, making sure that the ends are not immersed in order to avoid any electrical problems with the connectors. • Follow the disinfectant manufacturer's recommendation for contact time. • Rinse by partially immersing the cable in sterile or filtered water (0.2 µm) for 1 minute. Similar to the disinfection process, make sure that the ends are not immersed.
DRYING	• Dry with a single use lint-free towel.
MAINTENANCE, INSPECTION, AND TESTING	• Visually verify that there is no soiling, stains or residue on the cables. • Visually check that the cable is not damaged.
STORAGE	Ensure that required storage conditions are met. See section "CONDITIONS FOR OPERATION, TRANSPORT, AND STORAGE".
MANUFACTURER CONTACT	See manufacturer information in this document.

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