
NISSHA

MEDICAL TECHNOLOGIES

SWAROPLATE®

SPLIT | NON-SPLIT | TWIN-SAFE® | OVALPLATE™

Neutral Electrodes

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Swedish (SV)..... 16

Vietnamese (VI) 17

Instructions for Use Neutral Electrodes

This document applies to all SWAROPLATE neutral electrodes of NISSHA MEDICAL TECHNOLOGIES Ltd.

Intended Use

Neutral electrodes are an accessory for monopolar HF-surgery and represent the large area and low impedance contact with the patient's skin required for returning the electric current to the HF-generator. Neutral electrodes are self-adhesive, ready-to-use disposable products.

The SWAROPLATE neutral electrodes of Nissha Medical Technologies are not intended for use in High Current Mode. These are applications with high current and / or long activation periods where heating factors of more than 30 A2s in a 60 s period occur.

Indications

Indication is any medical condition that is treated by monopolar biterminal electrosurgery, except applications in high current mode.

Patient groups are neonates (<5kg), children (5-15kg) and adults (>15kg). Electrodes are specified for one or more of these patient groups.

Contraindications

Neutral electrodes must not be applied to injured or damaged skin. Therefore, a general contraindication is applicable to patients with large-area skin injuries (e.g. burns) if there is not enough intact skin application area available.

General Warnings Safety and Precautions

Neutral electrodes may be used by medical professionals only who must have undergone training regarding monopolar HF-surgery. Ignoring these instructions for use may cause severe damage to the patient's skin!

- The SWAROPLATE neutral electrodes of Nissha Medical Technologies must not be used during applications in High Current Mode! Because of the high HF-current and / or during the long activation of the HF-current in this mode the temperature under the neutral electrode increases more than at conventional procedures and this could lead to undesired tissue damage.
 - Keep activation times as short as possible. If long activation is necessary, allow time between activations to allow the tissue under the neutral electrode to cool.
 - The neutral electrodes are ready-for-use. Do not use additional contact gel.
 - Prior to application, make sure that the neutral electrode is intact and that the expiration date has not been exceeded.
 - Any modification of the neutral electrode (e.g. cutting) voids any liability of the manufacturer.
 - Neutral electrodes are single use products. Do not reapply an improperly applied electrode.
 - Review the application scopes for the different types of electrodes according to the weight and power specifications of the HF-generator as indicated on the packaging (max. 250W for children, max. 150W for newborns).
 - If the HF-generator has a contact quality monitoring system, (CQM, e.g. "REM," "NESSY," "EASY") always use split electrodes. Only when a split electrode is used with a CQM system, the contact with the patient is monitored automatically.
 - With non-split neutral electrodes, regularly check that the entire surface area of the neutral electrode has contact with the skin (especially after moving the patient) because this is not monitored by the HF-generator.
 - If the desired surgical effect is not achieved, stop the procedure and verify a good neutral electrode contact to the skin before proceeding with activation.
 - Use neutral electrodes that are compatible with your HF-generator only. Therefore, review the current Declaration of Compatibility (see Website).
 - Always observe the instructions for use of the HF-generator.
 - Any serious incident that has occurred in relation to the neutral electrodes should be reported to the manufacturer and the competent authority of the state in which the user and/or patient is established.
1. **Positioning** - Select an application site for the neutral electrode as close as possible to the operating site, but not closer than 20 cm. Avoid proximity to ECG electrodes and ensure that the heart is never in the current path. Preferred application sites are the thigh and upper arm. Always ensure the patient is not lying on top of the neutral electrode.
 2. **Avoid** - The neutral electrode must not be put on scars, inflamed or injured skin, bony parts of the body, above metal implants or above massive subcutaneous fat tissue.
 3. **Prep** - Carefully remove the patient's hair at the application site. If necessary, clean it. The application site must be dry and free of skin creams. Allow any used cleaning products to dry out completely.
 4. **Open** - Open the pouch and remove the neutral electrode. Once the pouch has been opened, use within 7 days. Always keep unused neutral electrodes in the pouch and close it by folding over.
 5. **Peel** - Check that the electrode and connected cable (if applicable) are free from any damages. Then, peel off the neutral electrode from the protective foil. Do not use damaged electrodes.
 6. **Caution** - Rectangular neutral electrodes must be applied in such a way that the long edge faces the site of surgery. Electrodes with an oval shape can be oriented in any direction due to their design.
 7. **Apply** - Attach the neutral electrode with its full area and without creases. Make sure that there are no air bubbles below the neutral electrode and that the protective edge (foam edge) is applied to the patient completely and without creases.
 8. **Connect** - Link the neutral electrode with the generator according to its instructions for use. Avoid contact of the cable with the patient and other cables.
 9. **Remove** - Remove the neutral electrode carefully from one corner after use to avoid skin irritations. Please note that disposal must be in line with the applicable laws and requirements of the respective country.

Explanation of Packaging Symbols

A. Electrode suited for adults (weight > 15 kg)	K. Temperature limitation	U. Quantity
B. Electrode suited for children (weight 5 kg to 15 kg)	L. Keep dry	V. Lot Number
C. Electrode suited for newborns (weight < 5 kg)	M. Keep Away from sunlight	W. Catalog Number
D. Electrode with integrated 3 m cable and REM connector	N. Green Dot scheme	X. Manufacturer
E. Electrode with integrated 3 m cable and NON-REM connector	O. CE Marking of Conformity and number of the involved Notified Body	Y. Authorized representative in the European Community
F. Electrode with integrated 3 m cable and 6.3mm jack connector	P. Medical Device	Z. Authorized Representative in Switzerland
G. For professional users only	Q. Date of Manufacture	AA. Importer
H. Consult instructions for use	R. Use by date	
I. Single use. Do not re-use	S. Does not contain PVC	
J. Not made with natural rubber latex	T. Do not use if package is damaged	

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Nissha Medical Technologies Ltd
Nissha Medical Technologies SAS