

The 4020/03 Pleated HEPA filter protects the patients airways effectively from exogenous microbial loads, reducing the risk of extrinsic colonization and infection.

FEATURES



LINE OF DEFENSE

Reduce the risk of extrinsic colonization and infection.



HYDROPHOBIC FILTER MEMBRANE

Cannot pass potentially contaminated fluids under normal conditions.



INHIBIT MICROORGANISMS

High bacterial and viral filtration efficiency rates to reduce risk of cross contamination.





UNIVERSAL SIZING

22mm x 22mm connector ports.



Bacterial and Viral Filtration to reduce Cross Contamination

FILTER PLACEMENTS IN A VENTILATOR SYSTEM



Filter can be used by itself or in line with other filters.

PACKAGING

POUCH

1 non-sterile filter per bag

CASE

50 Bags / case

4020/03 PLEATED HEPA FILTER

Part #	Filter Type	Dimensions	Efficiency	Connection
32031705	PLEATED HEPA FILTER	81.5mm x 68.5mm	BFE 99.9999% VFE 99.9999%	22M-22F