


Technical Documentation	Reusable Pressure Cables DECLARATION OF COMPATIBILITY NMT	GRD12_02-05-001
		REV: AA
		PAGE: 1 OF 2

REF		
64330	Philips	BD Argon
64474	Fukuda 8000 Series	BD Argon
64551	GE Marquette	BD Argon
64760	Draeger Siemens	BD Argon
82474	Fukuda 8000 Series	Utah Medical
82551	GE Marquette	Utah Medical
83330	Philips	PVB Codan
83474	Fukuda 8000 Series	PVB Codan
83551	GE Marquette	PVB Codan
83792	Spacelabs	PVB Codan
85330	Philips	B.Braun
85474	Fukuda 8000 Series	B.Braun
85551	GE Marquette	B.Braun
85761	Dräger	B.Braun
86112	Schiller CM8	Abbott
86330	Philips	Abbott
86474	Fukuda 8000 Series	Abbott
86481	SAMS	Abbott
86540	Mennen	Abbott
86551	GE Marquette	Abbott
86767	Draeger Siemens	Abbott
86792	Spacelabs	Abbott
87330	Philips	Smith Medex
87474	Fukuda 8000 Series	Smith Medex
87551	GE Marquette	Smith Medex
87560	Mindray	Smith Medex
87767	Draeger Siemens	Smith Medex
IP-88120	GE Critikon Dinamap	Edwards / BD
IP-88330B	Philips	Edwards / BD
IP-88330B-10M	Philips	Edwards / BD
IP-88330B-2M	Philips	Edwards / BD
IP-88330B-5M	Philips	Edwards / BD
IP-88330-V3B	Philips	Edwards / BD
IP-88470	Kontron	Edwards / BD
IP-88474	Fukuda 8000 Series	Edwards / BD
IP-88545	Nihon Kohden	Edwards / BD
IP-88550	GE Eagle Tram	Edwards / BD
IP-88551	GE Marquette	Edwards / BD
IP-88552	Bard	Edwards / BD
IP-88555	Abbott CentriMag	Edwards / BD

REF		
IP-88560	Mindray	Edwards / BD
IP-88700	Medtronic Physiocontrol	Edwards / BD
IP-88760	Draeger Siemens	Edwards / BD
IP-88780	GE Datex Cardiocap	Edwards / BD
IP-88792	Spacelabs	Edwards / BD

DOCUMENT HISTORY

REVISION	DATE	EDITOR	DESCRIPTION OF CHANGE
AA	2026-08-10	Daren Davies	Initial Release per DDP-0059, based on TESPI290W70010 Iss A