



Protect eyes and face against liquids and dust with LATITUDE™ Protective Face Shield. Manufactured in the USA with 3M material, our full face shield offers maximum protection on all sides of the face, protecting the eyes, nose and mouth from contact with liquids and dust. LATITUDE™ is not impact resistant nor intended for use against open flames or sparks.

Together we are addressing not only the needs of the medical community on the front-line fighting COVID-19 but also the needs of industry to support their reopening during the current pandemic.



LIGHTWEIGHT & OPTICALLY CLEAR

Provides clear, unobstructed view. Lightweight and clear. Effective for use against liquids and dust. 3M Shield Technology includes Anti-Fog coating.



MAXIMUM PROTECTION

Complete coverage including the sides of the face preventing all contact of liquids and dust from eyes, mouth, and nose.



COMFORTABLE AND FUNCTIONAL

LATITUDE™ comes with a comfortable black foam. It also easily goes over glasses or goggles. Latex, Fiberglass and Distortion Free.

LATITUDE - PROTECTIVE FACE SHIELD

Description

LATITUDE™ PLUS GENERIC HEADBAND
3M™ Anti-Fog Face Shield Film, 9960FS, 6.8 mil

Anti Fog



Per Bag

1

Per Case

200

MOQ

200

LATITUDE™ PLUS CUSTOM HEADBAND
3M™ Anti-Fog Face Shield Film, 9960FS, 6.8 mil
Custom 4-color image with logo



1

200

4800

HEADBAND OPTIONS

LATITUDE™ PLUS:

Black Label - Black Foam/Elastic



LATITUDE™ PLUS Custom:


4 Color Custom Artwork



LATITUDE™ PLUS - CUSTOM HEADBANDS

Increase brand awareness while keeping your staff and patrons safe.
Design Guidelines: <http://design.graphiccontrols.com>

1 Shield / Bag
200 Shields / Case



3-13-01-074 (4/2020)

