

Flexstone Basket

Single-Use Steerable Stone Retrieval Basket

FEATURES



REVOLUTIONARY STEERABILITY

The only basket offering independent steerability for precise navigation in challenging anatomies and difficult-to-reach stones.



BUILT FOR THE TOUGHEST CASES

Engineered to handle complex anatomies and lower pole stones while enabling faster, more efficient procedures in routine cases.



SEAMLESS INTEGRATION



Fully compatible with third-party scopes, ensuring effortless adoption into any clinical workflow.




PROTECT YOUR EQUIPMENT


Independent steerability reduces the steering burden on scopes, improving their longevity and minimizing wear and tear.

- World's First Steerable Kidney Basket
- Reduced Procedure Time
- Optimized Retrieval Efficiency






STONE CAPTURE TIME



Model	FS-A-115-21
Outer Diameter	0.7mm (2.1 Fr)
Basket Diameter	10 mm
Working Length	115 cm
Tip Steering	±90°



UNMATCHED CONTROL

Advanced zero-tip retrieval cage, 2.1F (0.7mm) polyimide transmission, ergonomic handle, and independent steerability ensure safe and effective stone extraction.

