

Pediatric Defibrillator Compatibility

Match Defibrillator Brands to Nissha Medical Part Numbers



BRAND / MODELS OF DEFIBRILLATOR	COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS
 Physio-Control LifePak - 15 LifePak - 20e	Pediatric-Physio
 Zoll Medical R-Series X-Series Propaq MD	Pediatric-Zoll

Pediatric Defibrillation Electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Please contact Nissha Medical Technologies at 1.800.669.6905 or HS.Info@NisshaMedical.com for any questions regarding electrode compatibility.