



## High Flow & Standard Flow Nasal Cannulas

### FEATURES



#### **DIVIDED CANNULA**

High flow nasal cannula used in conjunction with bubble humidifier delivers up to 15LPM



#### **PATIENT COMFORT**

Soft tubing reduces device-related pressure ulcers



#### **MULTIPLE SIZES AND APPLICATIONS**


Premature to pediatric sizing allows for system-wide standardization



#### **ECONOMICAL SOLUTION**

Direct economic benefits with reduction of wound care expenses, and elimination of need for foam pads and wraps

Reduce the risk of vascular soft tissue damage with the softness and elasticity of the cannula material.



6-Channel kink resistant tubing

Elastomeric soft tubing

#### **THREADED NUT CONNECTOR**

Universal connector eliminates the need for a christmas tree connector



The unique elastomeric properties of the over-the-ear Comfort Soft Plus® tubing minimize friction and shear forces on the skin at all contact points and act as a "strain relief" when a patient turns their head. As a result, Comfort Soft Plus cannulas are more comfortable, enhancing patient compliance and therapy outcomes.

## High Flow Heated Humidification Nasal Cannulas

### WITH 15MM OD AND 22MM ADAPTER, SOFT EARS EARPIECE

Part #	Ref	Description	Size	Flow Range	Color	QTY
700000952	0557	Comfort Soft Plus® Heated Humidification Nasal Cannula	Premature	Up to 7LPM	Yellow	10/case
700000953	0558	Comfort Soft Plus® Heated Humidification Nasal Cannula	Infant	Up to 7LPM	Light Blue	10/case
700000954	0559	Comfort Soft Plus® Heated Humidification Nasal Cannula	Neonatal	Up to 7LPM	White	10/case
700000955	0569	Comfort Soft Plus® Heated Humidification Nasal Cannula	Pediatric	Up to 8LPM	Green	10/case

## Comfort Soft Plus® Standard Flow Nasal Cannulas

### MINIMIZES EAR SORES, KINK RESISTANT TUBING, CONSISTENT OXYGEN DELIVERY

Part #	Ref	Description	Tubing Length	O2 Connector	Flow Range	QTY
<b>Premature Nasal Cannulas</b>						
700000943	0152	Comfort Soft Plus® Standard Flow Nasal Cannula	7'	Trumpet	Up to 2LPM	25/case
700000944	0153	Comfort Soft Plus® Standard Flow Nasal Cannula & Threaded Nut	7'	Thread Grip	Up to 2LPM	25/case
<b>Neonatal Nasal Cannulas</b>						
700000946	0163	Comfort Soft Plus® Standard Flow Nasal Cannula & Threaded Nut	7'	Thread Grip	Up to 3LPM	25/case
700000947	0167	Comfort Soft Plus® Standard Flow Nasal Cannula	7'	Trumpet	Up to 3LPM	50/case
<b>Infant Nasal Cannulas</b>						
700000948	0173	Comfort Soft Plus® Standard Flow Nasal Cannula & Threaded Nut	7'	Thread Grip	Up to 4LPM	50/case
700000949	0177	Comfort Soft Plus® Standard Flow Nasal Cannula	7'	Trumpet	Up to 4LPM	50/case
<b>Pediatric Nasal Cannulas</b>						
700000950	0183	Comfort Soft Plus® Standard Flow Nasal Cannula & Threaded Nut	7'	Thread Grip	Up to 6LPM	50/case
700000951	0187	Comfort Soft Plus® Standard Flow Nasal Cannula	7'	Trumpet	Up to 6LPM	50/case