

Defibrillation Electrodes

Defib Pads for all departments of the hospital



36 MONTH SHELF LIFE LONGEST ON THE MARKET

Extended shelf life reduces waste and saves money



60 INCH LEAD WIRES LONGEST ON THE MARKET

Increased flexibility with equipment placement & clinician mobility



ADVANCED HYDROGEL

Unique gel provides maximum adhesion



SMALLER PAD SIZE

Allows for easier placement without crowding



LARGER CONDUCTOR PLATE

Larger conductive area disperses energy more effectively



COLOR CODED PACKAGING

The correct product, for your device.

VISIBLE RELEASE LINER

Prevents slip hazards for staff.

INSTRUCTIONS FOR USE

User checklist, on the release liner.

COLOR CODED PADS

Anterior electrode is blue foam and the posterior electrode is white foam. This helps train staff to look for colors to aid in pad placement.

LEADS OUT OPTION

Allows pads to be pre-connected prior to opening pouch.

PHYSIO CONTROL



ZOLL





CARDIAC SCIENCE





Defibrillation Electrodes

Defib Pads for all departments of the hospital





PHYSIO CONTROL

Product #	Ref #	For Defibrillator Models	Leads	Pad & Leadwires
6600017H	C100AC-PHYSIO	Lifepak 15, 20/20e, 1000	In	
6600022H	C100LOAC-PHYSIO	Lifepak 15, 20/20e, 1000	Out	
6600093H	T100AC-PHYSIO	Lifepak 15, 20/20e, 1000	In	
6600100H	T100LOAC-PHYSIO	Lifepak 15, 20/20e, 1000	Out	



ZOLL MEDICAL

Product #	Ref #	For Defibrillator Models	Leads	Pad & Leadwires
6600019H	C100AC-ZOLL	R-Series, X-Series	In	
6600023H	C100LOAC-ZOLL	R-Series, X-Series	Out	
6600095H	T100AC-ZOLL	R-Series, X-Series	In	
6600102H	T100LOAC-ZOLL	R-Series, X-Series	Out	



CARDIAC SCIENCE

Product #	Ref #	For Defibrillator Models	Leads	Pad & Leadwires
6600088H	T100-CS	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	Can be used in place of Cardiac Science Part #9131-001



PEDIATRIC DEFIBRILLATION ELECTRODES




PHYSIO CONTROL

Product #	Ref #	For Defibrillator Models	Leads	Pad & Leadwires
6600068H	*PEDIATRIC-PHYSIO	Lifepak 15, 20e	In	



ZOLL MEDICAL

Product #	Ref #	For Defibrillator Models	Leads	Pad & Leadwires
6600070H	*PEDIATRIC-ZOLL	R-Series, X-Series, Propaq MD	In	



*HeartSync pediatric defibrillation electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

For additional information, defibrillator compatibility and usage visit: [HS.NisshaMedical.com/defib](https://www.hs-nissha.com/defib)