



## EMS Electrodes

12 Lead ECG Kit

### FEATURES



**V1-V2**

Small rectangular electrode with two sensors



**V3-V6:**

Small rectangular shaped electrode for placement in tight areas



**RA, LA, RR, LL:**


Large round electrodes for superior adhesion

Our single use pouch contains 10 electrodes perfect for EMS diagnostic stress tests.

### PACKAGING


**HYGIENIC, DISPOSABLE  
SINGLE-USE PACKS**

Reduces risk of cross contamination of unused electrodes



**CASE**

FDA / GS1 Compliant Labels




Built with time savings in mind, our single-use pouch contains 10 electrodes perfect for EMS diagnostic stress tests.

**REF #A10052-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT**  
**(PRODUCT # 6502771V) 10 ELECTRODES/ KIT | 50 KITS / CASE**

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.875" x 1.25" / 47.625mm x 31.75mm	Rectangle	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

## WHAT'S IN A KIT?



**4** 2.25" / 57.15mm  
Foam | Wet Gel



**4** 1.875" x 1.25" /  
47.625mm x 31.75mm  
Foam | Wet Gel



**1** 1.75" x 1.625" /  
44.45mm x 41.275mm  
Perforated | Foam | Wet Gel



Perforation allows separation when desired