

PerformancePlus™ Electrodes

Stress Test Kit

FEATURES



V1-V2

Small rectangular electrode with two sensors



V3-V6:


Small teardrop shaped electrode for placement in tight areas



RA, LA, RR, LL:

Large round electrodes for superior adhesion

Unique single-use pouch containing 10 electrodes. For use in diagnostic stress testing.



PACKAGING

HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



CASE
FDA / GS1 Compliant Labels




Our unique single use pouch containing 10 electrodes is perfect for diagnostic stress testing. Electrodes come in varying sizes for ease of placement in designated positions.

REF #A10015-10S - STRESS TEST KIT - FOR 10 LEAD PLACEMENT
(PRODUCT # 6502923V) 10 ELECTRODES/ KIT | 50 KITS / CASE

Size	Shape	Per Pouch	Wet Gel
2.25" / 57.15mm	Round	4	
1.5" x 38.1mm	Teardrop	4	
1.75" x 1.625" / 44.43mm x 41.275mm	Square, dual element	1	

WHAT'S IN A KIT?



4 2.25" / 57.15mm
Foam | Wet Gel



4 1.5" / 38.1mm Teardrop
Foam | Wet Gel



1 1.75" x 1.625" /
44.45mm x 41.275mm
Perforated | Foam | Wet Gel



Perforation allows separation when desired