



FEATURES

- Conductive adhesive gel
- Patented traction liner prevents tabs from sliding off the liner - saving time and money
- Unique liner card ensures packs do not brick together for ease of removal and application
- Pliable substrate conforms to the body
- Adhesive solid gel requires minimal clean-up
- Latex-free

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

Latex-free Solid Gel



SilveRest™ Tab Electrodes

Diagnostic Resting EKG Tab Electrode

Vermed SilveRest™ Tab Electrodes for resting EKG provide excellent adhesion and quality traces. All of our tab electrodes meet or exceed the AAMI standard for disposable pre-gelled EKG electrodes.


A10022 - .984" X .90" / 24.99MM X 22.86MM - AG/AGCL

- Specially designed strain relief tab is pliable, allowing electrode to remain secure.
- Ultra-low impedance (<100 ohms) promotes rapid baseline stabilization.

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10022	10	100	1000	4000	
A10022-1000BX	10	100	n/a	1000	

A10201 - 1.34" X .90" / 34.036MM X 22.86MM - C/AGCL

- Proprietary low residue gel for easy clean-up.
- Larger clip attachment area uses small or large clips without compromise.

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10201	10	100	1000	4000	
A10201-500	10	100	500	3000	
A10201-1000BX	10	100	n/a	1000	

A10023 - 1.125" X 1.438" / 28.575MM X 36.525MM - AG/AGCL

- Proprietary adhesive solid gel for optimal adhesion.
- Larger surface area contact which reduces impedance and provides superior adhesion.
- Repositionable and packaged in resealable pouches

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10023	10	100	1000	4000	
A10023-1000BX	10	100	n/a	1000	

A10022
Tab | Solid Gel
Silver / Silver



A10201
Tab | Solid Gel
Silver / Carbon



A10023
Regal Tab | Solid Gel
Silver / Silver

