


SilveRest™ 10-Pack Tab Electrodes

Diagnostic Resting EKG Tab Electrode



FEATURES



TRACEABILITY

Each single-use pack is lot coded and has expiration dates.



HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Eliminates risk of cross contamination of unused electrodes.



IMPROVED FUNCTIONALITY

Eliminates electrode failure due to improper sealing and storage of unused electrodes.



MR CONDITIONAL TAB ELECTRODES


Can be worn (but not used) during MRI procedures.

Introducing the market's first single use 10-pack of tab electrodes.

PACKAGING

HYGIENIC, DISPOSABLE SINGLE-USE PACKS

Reduces risk of cross contamination of unused electrodes



INNER BOX

Single-use packs are packaged in inner boxes within the case

CASE

FDA / GS1 Compliant Labels



SilveRest™ 10 Pack Tab Electrodes

Diagnostic EKG Resting Electrode

The Joint Commission (TJC) is mandating multi-use packages be documented for first use and expiration dates. Single-use packs eliminate the need to document these dates. Single-use pack electrodes also provide a solution to assist in compliance with IC.02.02.01 standard for infection control.

A10202-10 - 34mm x 23mm - C/AgCL



Product #	Ref #	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6503179V	A10202-10	10	500	4000		 < 24 Hours