



FOR USE IN NICU, DIAGNOSTICS, APNEA, AND PICU



FEATURES

- **Infant Breathable Electrodes:** Designed for patient comfort. Gentle adhesive and low chloride hydrogel provide a conforming and non-irritating electrode. Excellent for use in APNEA monitoring and all other pediatric applications.
- **Soft Cloth Electrodes:** Adhesive solid gel electrodes with soft cloth material are available in various sizes to suit your needs.
- **Foam Pediatric Electrodes:** Our foam pediatric electrodes are available in 1" diameter to allow for easy placement and features solid gel that is gentle on the skin.

1.800.669.6905 | Vermed.com

Electrodes • Chart Paper • Defib Pads • Cables & Leadwires

Latex-free



Wet Gel



Solid Gel




TenderTrode™


Tender Care for Tender Patients

Vermed's TenderTrode™ line represents a complete selection of electrodes for infant and pediatric use. Each electrode is carefully designed with gentle, hypoallergenic adhesives and soft conforming materials to ensure comfort and effectiveness on young patients.


A10001 - 1.625" / 41.275MM - ROUND TAPE ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Wet Gel
A10001	3	3	n/a	600	
A10000-60	3	60	n/a	600	


A10040 - 0.875" / 22.225MM - SQUARE CLOTH ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10040-60	3	60	n/a	600	
A10041-5	5	5	n/a	600	
A10041-6	6	6	n/a	600	

A10041 - 0.875" X 1.375" / 22.225MM X 34.925MM - RECT. CLOTH ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10041-60	3	60	n/a	600	
A10041-5	5	5	n/a	600	
A10041-6	6	6	n/a	600	

A10042 - 1.375" / 34.925MM - SQUARE CLOTH ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10042-60	3	60	n/a	600	
A10042-10	5	10	n/a	450	

A10050 - 1" / 25.4MM - ROUND FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10050-60	3	60	n/a	600	

P7-3 - 0.875" / 22.225MM - TEAR DROP FOAM ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10050-60	3	3	108	1080	


A10027 - 0.5" X 0.9 / 12.7MM X 22.86MM - AG/AGCL TINYTAB ELECTRODE

Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel
A10027	10	100	n/a	1000	

A10001
Tape | Wet Gel



A10040
Cloth | Solid Gel



A10041
Cloth | Solid Gel



A10042
Cloth | Solid Gel



A10050
Foam | Solid Gel



P7-3
Foam | Solid Gel



A10027
Tab | Solid Gel
Silver / Silver

