



NMT Ultra Thermal Films

Nissha Medical Technologies (NMT) Ultrasound films are designed to be utilized in medical diagnostic applications with use of standard ultrasound systems. When comparing your current film to NMT Ultra you will find the physical attributes to be the same as other leading manufacturers.

In our recent study comparing NMT Ultra Film to other brands, we found that NMT not only exceeds expectations as a cost-effective alternative, but also stands out in key performance criteria such as graininess, contrast, and density. We believe NMT Ultra delivers exceptional value through its high quality and low cost per print.

HIGH QUALITY, COST EFFECTIVE ULTRASOUND FILM ALTERNATIVE


Dynamic Thermal Response



Graininess Evaluation



SONY UPP-110HD



NMT ULTRA 110HD

32019066 - NMT ULTRA 110HG THERMAL VIDEO PAPER

Product #	Per Box	Per Case	UOM
32019066	5	50	Roll

32019065 - NMT ULTRA 110HD THERMAL VIDEO PAPER

Product #	Per Box	Per Case	UOM
32019065	5	50	Roll

ULTRASOUND FILMS

NISSHA
MEDICAL TECHNOLOGIES





NMT THERMAL FILM



OTHER COMPETITOR

- ✓ High Humidity and Heat Resistance for long storage
- ✓ Top Coat Layer optimally matches the printer heads of printer for continuous printing
- ✓ Easy Tearing along the cross direction
- ✓ Water Resistant
- ✓ Anti-Static Layer
- ✓ Minimal Curling
- ✓ High Quality Prints



NMT ULTRA 110

OTHER COMPETITOR

Roll Dimensions	110 mm x 18 m	110 mm x 18 m
Size	A6	A6
Packaging	5 Rolls / Box 50 Rolls / Case	10 Rolls / Box
Feature	High Gloss	High Gloss
Print	Black & White	Black & White
Applications / Usage	Mainly Ultrasound Applications	Mainly Ultrasound Applications
Sony Printer Compatibility	UP-X898MD / D898MD / D898DC / 897MD / D897	UP-X898MD / D898MD / D898DC / 897MD / D897
Mitsubishi Printer Compatibility	P95 / P93	N/A
Storage Conditions	5°C/41°F ~ 30°C/86°F	5°C/41°F ~ 30°C/86°F

Roll Dimensions	110 mm x 20 m	110 mm x 20 m
Size	A6	A6
Packaging	5 Rolls / Box 50 Rolls / Case	10 Rolls / Box
Feature	High Density	High Density
Print	Black & White	Black & White
Applications / Usage	Ultrasound, Dental, and Microscopy Applications	Ultrasound, Dental, and Microscopy Applications
Sony Printer Compatibility	UP-X898MD / D898MD / D898DC / 897MD / D897	UP-X898MD / D898MD / D898DC / 897MD / D897
Mitsubishi Printer Compatibility	P95 / P93	N/A
Storage Conditions	5°C/41°F ~ 30°C/86°F	5°C/41°F ~ 30°C/86°F

SIDE BY SIDE COMPARISON

NMT ULTRA 110HD



Competitor 1 10HD



NMT ULTRA 110HG



Competitor 1 10HG



NISSHA

MEDICAL TECHNOLOGIES

3-13-01-083-BUF AB (09/2025)

NORTH AMERICA • BUFFALO, NY | SANDPOINT, ID | COLORADO SPRINGS, CO | NEW RICHMOND, WI | NORWALK, OH | GANANOQUE, ONTARIO
CARRIBEAN • DOMINICAN REPUBLIC
EUROPE • TORBAY, UNITED KINGDOM | DUREN, GERMANY | RUM, AUSTRIA | CERGY, FRANCE