

Pediatric Defibrillator Compatibility

Match Defibrillator Brands to Nissha Medical Part Numbers



BRAND / MODELS OF DEFIBRILLATOR		COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS
	Philips	
	HeartStart XLT	
	HeartStart XL+	
	HeartStart XL	
	MRX	
	HeartStart Intrepid	Pediatric-Philips
	Physio-Control	
	LifePak - 12	
	LifePak - 15	
	LifePak - 20 / 20e	Pediatric-Physio
	Zoll Medical	
	R-Series	
	X-Series	
	M-Series	
	E-Series	
	Propaq MD	Pediatric-Zoll

Pediatric Defibrillation Electrodes are NOT compatible with AED's. AED's cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Please contact Nissha Medical Technologies at 1.800.669.6905 or HS.Info@NisshaMedical.com for any questions regarding electrode compatibility.