


PATIENT MONITORING



EUROPE

NISSHa
MEDICAL TECHNOLOGIES

TABLE OF CONTENTS



ECG ELECTRODES

Product Packaging	4	Prewired Paediatric	16
Electrode Anatomy	6	SilveRest Tabs	20
VersaTrode™	8		
Swaromed	12		
Swaromed Paediatric	14		



DEFIB PADS

Defib Pad Anatomy	22	Zoll Medical	24
Philips Medical	24	Cardiac Science	24
Physio-Control	24		



GROUNDING PADS (Diathermy Pads | Neutral Electrodes)

Swaroplate™ Split Pads	27	Swaroplate™ Non-Split Pads	27
Swaroplate™ Split - Prewired	27	Swaroplate™ Non-Split - Prewired	27


Let Colour be Your Guide

We know that selecting the best product for every application saves you time and money - so we make it easy. Colour coding for quick reference means confidence in your selection. Every section through out the catalog has its own unique colour, allowing you to quickly and easily find exactly what you are looking for.




GROUNDING PADS (Continued)

TwinSafe™ Grounding Pads	28	OvalPlate™ Non-Split - Prewired	28
TwinSafe™ - Prewired	28	Grounding Pad Anatomy	29
OvalPlate™ Split Pads	28	Grounding Pad Cables	29
OvalPlate™ Split - Prewired	28	Grounding Pad Packaging	29
OvalPlate™ Non-Split Pads	28		



MEDICAL CHART PAPERS

Burdick	31	Instrumentation Labs	33	Quinton	34
Corometrics	31	Marquette GE	33	Radiometer	34
Datascope	31	Medtronic	33	Schiller	34
Dinamap	31	Mennen	33	Siemens	35
Edan	32	Micro Medical	33	Sonicaid (Huntleigh)	35
Electrodyne	32	Mindray	33	Spacelabs	35
Fukuda Denshi	32	Miscellaneous	33	St. Johns	35
GE	32	Mortara	34	Steris	35
Getinge	32	Nihon Kohden	34	Welch Allyn	35
Hellige	32	Philips	34	Vitalograph	35
Hewlett Packard	32	Physio-Control	34	Zoll	35
Humphrey	33	Protocol	34		




FACILITIES CHARTS

Helmer	37	United Electric	37
Jewitt	37	Jules Richard	37



ULTRASOUND GELS

Vermed	39
Aquasonic 100	39
Aquasonic Clear	39
SCAN Gels	40
Ultrasound Probe Covers	40
Thermasonic Gel Warmers	40
Conductive Media	40



FETAL BELTS

Abdominal Transducer Belts	41
----------------------------	----



ULTRASOUND FILM

Sony	43
Mitsubishi	44
Vermed Ultra	45

ELECTRODE PACKAGING

BULK POUCHES

25 or more electrodes.

BULK POUCHES (Front)

Country Specific Languages

BULK POUCHES (Back)

Pictogram and Directions
For Use (DFU's) in Country
Specific Languages




CASE LABELS

Product # and information
placed on case packaging

INDIVIDUAL ROLLS

Individually shrink wrapped rolls bulk packed within the case




MULTIROLL PACK

Sets of 5 or 10 rolls shrink wrapped together




INNER SLEEVES

5 individual rolls are packaged in inner sleeves within the case



FAN FOLD

Individually shrink wrapped fanfold packs within a case



ULTRASOUND FILM CASE

Individually packaged rolls of film packed within the case





STAINLESS SNAP

Stainless steel snap stops corrosion from occurring, eliminating artifact in the tracing. Contains a Silver / Silver Chloride sensor.

SUBSTRATE


Vermed Electrodes are one millimeter thick. Our foam electrodes provide more strength and less stretch during placement.

WET GEL

Our wet gel electrodes have been developed to ensure a perfect chloride content for more reliable traces.

LINERS

Vermed wet gel electrodes come on a clear liner that has a unique protective gel cap that protects the integrity of the gel and prevents its dispersion on the adhesive area



GEL OPTIONS



Solid Gel



Wet Gel

LABEL

The dot label makes sure the snap is held in place on the electrode. It is also generally colour coordinated to help identify product lines and application.

CARBON SNAP

Our Carbon snaps are MR Conditional. It can be worn during MRI Scans, but not actively in use. Radiolucent and contains a Silver / Silver Chloride sensor.

ADHESIVE

Our adhesive material is made from a medical grade acrylic - ensuring that your electrodes stay in place. Larger electrodes allow for an increase in adhesion.

SOLID GEL

Solid gel is filled in at the depth of the foam to allow for smooth, flat application. Solid gel aids in the adhesion of the electrode and does not require clean up after use.



ALL MATERIALS



SNAP OPTIONS



Stainless Steel



Carbon Filled

MR Conditional



SUBSTRATE OPTIONS

Foam Electrodes

1mm Polyethylene Foam

Cloth Electrodes

Spunlace Nonwoven Fabric



VersaTrode™



Vermed's VersaTrode™ line of disposable ECG electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications. Each VersaTrode™ electrode is designed to work in all procedures: Stress Test, Long-Term Monitoring, Short Term Monitoring or Resting ECG.

FEATURES



WET & SOLID GEL ELECTRODES

Wet gel provides quick reading and stabilization. Solid gel improves adhesion and requires minimal clean up.



HYPOALLERGENIC ADHESIVE

Hypoallergenic adhesive and not made with natural rubber latex.



HIGH TACK OPTION

Offers stronger adhesive

VersaTrode™

General Monitoring Electrodes



60mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000937N	FW60	25	250	3000	Stainless Steel	Regular	



50mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000936N	FW50	50	400	4800	Stainless Steel	Regular	
6000935N	FW50R	50	400	4800	PP Carbon Blend	Regular	
6000941N	FS50	50	400	4800	Stainless Steel	Regular	
6000950N	FS50R	50	400	4800	PP Carbon Blend	Regular	
6001492N	HTFW50R	50	400	4800	PP Carbon Blend	High Tack	
6001493N	HTFW50	50	400	4800	Stainless Steel	High Tack	
6001491N	HTFS50	50	400	4800	Stainless Steel	High Tack	



36 x 50mm - Oval - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000972N	HTFW3654	25	250	3000	4mm	High Tack	
6000982N	HTFW3654	50	400	4800	4mm	High Tack	
6000940N	HTFS365	50	500	6000	Stainless Steel	High Tack	
6000933N	HTFW365	50	500	6000	Stainless Steel	High Tack	
6000982N	HTFW365B	30	300	3600	Adaptor	High Tack	



32mm - Square - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000930N	HTFW32	50	500	6000	Stainless Steel	High Tack	
6001023N	HTFS323	50	500	6000	Stainless Steel	Regular	
6001126N	FS32P	50	500	6000	PP Carbon Blend	Regular	
6001494N	HTFW323R	50	400	4800	PP Carbon Blend	High Tack	
6000930N	HTFW323	50	500	6000	Stainless Steel	High Tack	
6001200N	HTFW323B	50	400	4800	Adaptor	High Tack	



32mm - Teardrop - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000931N	FW32P	50	500	6000	Stainless Steel	Regular	
6001409N	4051	50	500	6000	Stainless Steel	High Tack	

Banana
Plug



Latex-free



Wet Gel



Solid Gel MR Conditional





Swaromed ECG Electrodes are pre-gelled (single-use) and have a Silver/Silver Chloride sensor. All adhesive foams used are equipped with a hypoallergenic medical-grade adhesive that is non-irritating to the skin.

FEATURES



SILVER/SILVER CHLORIDE SENSOR

High quality Ag/AgCL sensor helps produce clear traces.



HYPOALLERGENIC ADHESIVE

Hypoallergenic adhesive and not made with natural rubber latex.



STRICT QUALITY CONTROL DURING PRODUCTION

Complete confidence in quick, accurate traces

DEFIBRILLATION RESISTANT

XXXX - Using as a placeholder because is this actually true?



QUICK RECOVERY GEL

Wet gel provides quick reading and stabilization. Solid gel improves adhesion and requires minimal clean up.

Swaromed™

Adult General Purpose Electrodes



60mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000003C	1001	25	250	3000	Stainless Steel	Regular	●
6001122C	4001	25	250	3000	Stainless Steel	Regular	●
6000003N	1001	25	250	3000	Stainless Steel	High Tack	●



50mm - Round - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000018C	1009	50	400	4800	Stainless Steel	Regular	●
6000388N	1009A	25	250	3000	Stainless Steel	Regular	●
6000582C	4009	50	400	4800	Stainless Steel	Regular	●
6000436C	1010	50	400	4800	PP Carbon Blend	Regular	●
6000981C	4010	50	400	4800	PP Carbon Blend	Regular	●
6000665C	1012	50	400	4800	PP Carbon Blend	High Tack	●
6000341C	1019	50	400	4800	Stainless Steel	High Tack	●
6000616C	4019	50	400	4800	Stainless Steel	High Tack	●
6001481C	1085	25	250	3000	Adaptor	Standard	●



36 x 50mm - Oval - Foam Electrodes

Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000045C	1036	50	500	6000	Stainless Steel	High Tack	●
6000679N	1036A	50	500	6000	Stainless Steel	High Tack	●
6000578C	4036	50	500	6000	Stainless Steel	High Tack	●
6000398C	1040	50	500	6000	PP Carbon Blend	High Tack	●
6000069N	1084	25	250	3000	Stainless Steel	High Tack	●
6000609C	1084	30	300	3600	Stainless Steel	High Tack	●
6001237C	4040	50	500	6000	Stainless Steel	High Tack	●

Abrader Cap



Banana Plug



Latex-free



Wet Gel



Solid Gel



MR Conditional





Swaromed™ Paediatric

FEATURES

- Smaller electrodes for easier placement
- Available in wet or solid gel options
- Doodle on label helps smaller patients feel more comfortable
- Available with stainless steel or radiotranslucent PP Carbon snaps
- Banana plug adapters available

32mm - Teardrop - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000612C	1057	30	300	3600	Stainless Steel	Regular	
6001479C	4050	50	500	6000	Stainless Steel	Standard	
6001418C	4051	50	500	6000	Stainless Steel	High Tack	
6001279C	4057	25	250	3000	Stainless Steel	Standard	

ADULT & PAEDIATRIC 32mm - Square - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type
6000044C	1060	50	500	6000	PP Carbon Blend	High Tack	
6000588C	4060	50	500	6000	PP Carbon Blend	High Tack	
6000065C	1066	50	500	6000	Stainless Steel	High Tack	
6000534N	1066A	50	500	6000	Stainless Steel	High Tack	
6000581C	4066	50	500	6000	Stainless Steel	High Tack	
6000068N	1083	25	250	300	Stainless Steel	High Tack	
6000570C	1083	50	400	4800	Stainless Steel	High Tack	

Abrader Cap



Banana Plug



Latex-free



Wet Gel



Solid Gel



MR Conditional



Prewired Electrodes



FEATURES

NMT offers a complete line of prewired electrodes for use in the PICU/NICU and other neonatal applications. Available in various sizes including radiolucent leadwire options.




SOFT CLOTH ELECTRODES

Each electrode is carefully designed with gentle, hypoallergenic adhesives and soft conforming materials to ensure comfort and effectiveness on every young patient.

Prewired Electrodes

For Use in the PICU/NICU and Other Neonatal Applications

Prewired Neonatal ECG Electrodes



Product #	Description	Material	Per Bag	Per Box	Solid Gel	MR Conditional
40554IP	23 x 23 x 1mm Square	Cloth	3	150		
40556IP	23 x 23 x 1mm Square 1.5mm connector	Cloth	3	150		
40557IP	23 x 23 x 1mm Square 4mm connector	Cloth	3	150		
40558IP	23 x 23 x 1mm Oval	Cloth	3	150		

MR Conditional



Latex-free



Solid Gel



Tab Electrodes



Tab Electrodes by NMT for resting EKG provide excellent adhesion and quality traces. All of our tab electrodes meet or exceed the AAMI standard for disposable pre-gelled EKG electrodes.

FEATURES



MINIMAL CLEANUP

Conductive adhesive solid gel requires minimal cleanup.



COMFORTABLE AND FUNCTIONAL

Pliable substrate conforms to the body.



CONVENIENT PACKAGING

Patented traction liner prevents tabs from sliding off the liner and ensures packs do not brick together for ease of removal and application.



HYPOALLERGENIC

Tab electrodes are not made with natural rubber latex.

Tab Electrodes

Diagnostic Resting EKG Tab Electrode

A10022 - .984" x .90" / 24.99mm x 22.86mm - Ag/AgCL



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500043V	A10022	10	100	1000	4000		24 Hours

A10201 - 1.34" x .90" / 34.036mm x 22.86mm - C/AgCL



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500082V	A10201	10	100	1000	4000		24 Hours
6500083V	A10201-500	10	100	500	3000		24 Hours

A10023 - 1.125" x 1.438" / 28.575mm x 36.525mm - Ag/AgCL



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6500044V	A10023	10	100	1000	4000		24 Hours
6502902V	A10023-1000BX	10	100	n/a	1000		24 Hours
6502971V	A10023-10	10	10		4000		24 Hours

A10202-10 - 1.297" x .864" / 34mm x 23mm - C/AgCL



Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6503179V	A10202-10	10	10	500	4000		24 Hours


A10027 - 0.5" x 0.9 / 12.7mm x 22.86mm - AG/AGCL TinyTab



TinyTab electrodes are specially designed to fit on smaller pediatric patients

Part #	Product #	Per Strip	Per Pouch	Per Box	Per Case	Solid Gel	Duration of Use
6502930V	A10027	10	100	n/a	1000		24 Hours





60" / 1.5M LEAD WIRES

The longest on the market.
Provides more flexibility when
placing the defibrillator.



COLOR CODED PADS

RA electrode is white foam and the
LL electrode is blue foam. This helps
train staff to look for colors to aid in
pad placement.

SHELF LIFE

HeartSync Defib Pads have a 36 month shelf life
(3 years). That's approximately 50% longer than
our leading competitors.

POUCH OPENING

HeartSync Defib Pads have a traditional notched
tear opening, or an easy to peel alternative,
providing ease of use for clinicians.

DEFIB PACKAGING

Color coded pouches indicate defibrillator
equipment for easy identification. Leads out
option allows pads to be pre-connected
prior to opening pouch.

PEDIATRIC POUCHES

Light blue labels help eliminate user error
by providing fast and simple identification.

The OEM trademarks identified herein are the trademarks of the respective OEMs, and not NMT. NMT disclaims any affiliation, connection or association between its products and those of the respective OEMs.

PAD SIZE

Our pads are smaller than competitor pads allowing for easier placement on all body shapes and sizes.

CONDUCTOR PLATE

Our proprietary aggressive adhesive conductive hydrogel reduces the need for a large foam border allowing for a larger conductor plate than the competition.

INSTRUCTIONS FOR USE

User instructions and a checklist is on the release liner. Opaque liner does not become a slipping hazard for the staff on the treatment floor.

Instructions for use:

1. Remove electrodes from peel pouch. ☒
2. Ensure skin site is clean and dry. ☐
3. Remove electrodes from release liner. ☐
4. Apply electrodes firmly to patient's skin. ☐
5. Connect electrodes to defibrillator or patient cable. ☐
6. Connect ECG leads for demand pacing. ☐

PEEL
HERE
ONLY

PEEL
HERE
ONLY



Blue
Paediatric
Label

PHYSIO
CONTROL

PHILIPS
MEDICAL

CARDIAC
SCIENCE

ZOLL
MEDICAL

SCHILLER

NIHON-
KOHDEN

Defibrillation Electrodes



FEATURES



36 MONTH SHELF LIFE

LONGEST ON THE MARKET

Extended shelf life reduces waste and saves money



60 INCH LEAD WIRES

LONGEST ON THE MARKET

Increased flexibility with equipment placement & clinician mobility



ADVANCED HYDROGEL

Unique gel provides maximum adhesion



SMALLER PAD SIZE

Allows for easier placement without crowding




LARGER CONDUCTOR PLATE


Larger conductive area disperses energy more effectively

Defibrillation Electrodes

Defib Pads for All Departments of the Hospital




HeartSync's defib pads were designed with your needs in mind. Our defib pad line is not just different - it is a step above the other brands.




PHILIPS MEDICAL

Product #	Cross Ref #	For Defibrillator Models	Leads	Visibility
D311	C100AC-Philips	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
D309	T100AC-Philips	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL, HeartStart FR2+, HeartStart Forerunner	In	
D317	T100LOAC-Philips	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	Out	
D321	PEDIATRIC-Philips	XLT, MRX, HeartStart XL, MRX ALS, HP Codemaster XL	In	




PHYSIO CONTROL

Product #	Cross Ref #	For Defibrillator Models	Leads	Visibility
D304	C100AC-Physio	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
D313	T100AC-Physio	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	In	
D315	T100LOAC-Physio	Lifepak 9P, 15, 10C, 20, 12, 500, 500 DPS, 1000	Out	
D319	PEDIATRIC-Physio	Lifepak 9P, 15, 10C, 20, 12	In	



ZOLL MEDICAL

Product #	Cross Ref #	For Defibrillator Models	Leads	Visibility
D303	C100AC-Zoll	M-Series, E-Series, R-Series, X-Series	In	
D312	T100AC-Zoll	M-Series, E-Series, R-Series, X-Series	In	
D314	T100LOAC-Zoll	M-Series, E-Series, R-Series, X-Series	Out	
D318	PEDIATRIC-Zoll	M-Series, E-Series, R-Series, X-Series	In	



CARDIAC SCIENCE

Product #	Cross Ref #	For Defibrillator Models	Leads	Visibility
D326	T100-CS	PowerHeart AED G3 Plus, PowerHeart AED G3 Pro	Out	Can be used in place of Cardiac Science Part #9131-001

[†]Blue Highlighted part denotes paediatric.
For additional information, defibrillator compatibility and usage visit: HS.NisshaMedical.com/defib



Solid Gel



Grounding Pads



From our patented TwinSafe design to our standard SwaroPlate offering, NMT - HS Grounding (Diathermy) Pads are compatible with many electrosurgical generators. Pads are produced with the highest quality materials and designed to ensure complete contact during high-frequency surgery.

FEATURES



FUNCTIONALITY

Designed for Optimal Performance



EFFICIENCY

Enables Efficient Running of the Endoscopic Suite



INFECTION CONTROL

Minimises Infection Risk to Improve Patient Outcomes



COMPATIBILITY

Comprehensive Range to Simplify Procurement & Stock Management

Grounding Pads

Diathermy Pads/Neutral Electrodes

SWAROPLATE™ SPLIT PADS

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance @ 600kHz	Temp. Rise @ 60s
1211H	Adults	201 x 106mm	200 cm ²	144 cm ²	1 Ohm	< 3°C/700mA
1212H	Adults	206 x 122mm	220 cm ²	160 cm ²	1 Ohm	< 3°C/700mA
1213H	Adults	230 x 123mm	232 cm ²	169 cm ²	1 Ohm	< 3°C/700mA
1222H	Children	154 x 88mm	124 cm ²	71 cm ²	2 Ohm	< 5°C/500mA
1226H	Children	164 x 124mm	105 cm ²	71 cm ²	2 Ohm	< 5°C/500mA
1223H	Infants & Prematures	165 x 74mm	66 cm ²	32 cm ²	4 Ohm	< 5°C/350mA

SWAROPLATE™ SPLIT - PREWIRED


Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Base Pad	Connector
1251H	Adults	201 x 106mm	200 cm ²	144 cm ²	1211H	Valleylab REM
1255H	Adults	206 x 122mm	220 cm ²	160 cm ²	1212H	Valleylab REM
1252H	Children	154 x 88mm	124 cm ²	71 cm ²	1222H	Valleylab REM
1257H	Children	164 x 124mm	105 cm ²	71 cm ²	1226H	Valleylab REM
1253H	Infants & Prematures	165 x 74mm	66 cm ²	32 cm ²	1223H	Valleylab REM

SWAROPLATE™ NON-SPLIT PADS

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance @ 600kHz	Temp. Rise @ 60s
1201H	Adults	215 x 104mm	207 cm ²	134 cm ²	1 Ohm	< 3°C/700mA
1202H	Children	144 x 90mm	115 cm ²	72 cm ²	2 Ohm	< 5°C/500mA
1203H	Infants & Prematures	165 x 74mm	66 cm ²	35 cm ²	4 Ohm	< 5°C/350mA

SWAROPLATE™ NON-SPLIT - PREWIRED

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Base Pad	Connector
1261H	Adults	215 x 104mm	207 cm ²	134 cm ²	1201H	Valleylab NON-REM
1271H	Adults	215 x 104mm	207 cm ²	134 cm ²	1201H	Jack Connector



SWAROPLATE SPLIT
Grounding Pads



SWAROPLATE NON-SPLIT
Grounding Pads

Twinsafe™

TWINSAFE™ PADS

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Max Power Setting	Temp. Rise @ 60s
1280H	Adults & Children	176 x 122mm	169 cm ²	110 cm ²	n/a	< 5°C/350mA
1283H	Infants & Prematures	181 x 76mm	77 cm ²	37 cm ²	150 W	< 4°C/350mA

TWINSAFE™ - PREWIRED

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Max Power Setting	Connector
1285H	Adults & Children	176 x 122mm	169 cm ²	110 cm ²	N/A	Valleylab REM
1286H	Adults & Children	176 x 122mm	169 cm ²	110 cm ²	N/A	Jack Connector
1288H	Infants & Prematures	181 x 76mm	77 cm ²	37 cm ²	150 W	Valleylab REM

Ovalplate™

OVALPLATE™ SPLIT PADS

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance @ 600kHz	Temp. Rise @ 60s
1220H	Adults & Children	160 x 122mm	148 cm ²	105 cm ²	1.5 Ohm	< 4°C/700mA

OVALPLATE™ SPLIT - PREWIRED


Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Base Pad	Connector
1259H	Adults & Children	160 x 122mm	148 cm ²	105 cm ²	1220H	Valleylab NON-REM
1269H	Adults & Children	160 x 122mm	148 cm ²	105 cm ²	1220H	Jack Connector

OVALPLATE™ NON-SPLIT PADS


Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Impedance @ 600kHz	Temp. Rise @ 60s
1210H	Adults & Children	155 x 126mm	144 cm ²	105 cm ²	1.5 Ohm	< 4°C/700mA

OVALPLATE™ NON-SPLIT - PREWIRED

Cross Ref.	Use	Size (Length/Width)	Adhesive Surface	Conductive Surface	Base Pad	Connector
1261H	Adults & Children	155 x 126mm	144 cm ²	105 cm ²	1210H	Valleylab NON-REM
1271H	Adults & Children	155 x 126mm	144 cm ²	105 cm ²	1210H	Jack Connector



TWINSAFE
Grounding Pads



OVALPLATE SPLIT
Grounding Pads

OVALPLATE NON-SPLIT
Grounding Pads

Grounding Pad Anatomy

HYDROGEL TECHNOLOGY

Strong but gentle hydrogel adhesive provides immediate and steady impedance upon skin contact.

ADHESIVE FOAM BORDER

The highly flexible adhesive foam border and the conductive gel ensures that the electrode will stay secure.

REDUCED EDGE EFFECT


No heavily concentrated areas of current thus reducing the chance of patient burns.

SPLIT ELECTRODE

Compatible with most High Frequency (HF) Surgery Devices with Contact Quality Monitoring (CQM).

EXCELLENT SKIN COMPATIBILITY

Gel and adhesive foam are tolerated even by sensitive skin. There is no gel residue left on the skin after application.



Reusable Connection Cables

JACK CONNECTOR:

REF 1233: Length 3m
REF 1234: Length 5m
Connects to HF Surgery Device

VALLEYLAB REM:

REF 1239: Length 5m
Connects to HF Surgery Device

VALLEYLAB NON-REM:

REF 1236: Length 5m
Connects to HF Surgery Device

BLUE CLIP:

Connects to Grounding Pad

Packaging

POUCH

1 Electrode/Pouch



BOX

- 50 Neutral Electrode Pouches / Box
- 8 Prewired Neutral Electrode Pouches / Box





Medical Chart Paper

MEDICAL CHART PAPER

NMT is the ORIGINAL Chart Paper manufacturer, offering the widest assortment of chart papers for all departments in the hospital - from cardiology to the blood bank. We can consolidate the greatest number of suppliers for you – saving time and money. NMT manufactures in the UK and USA to strict OEM specifications. We offer the exact same grade of paper, dimensions and archival qualities as the Original Equipment Manufacturers. In fact, we make charts for many OEMs, offering the best **VALUE** product on the market - **OEM quality at a significant savings.**

Medical Chart Paper

Standardization Made Simple™

Burdick

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017246	BDK 007868	8.438" x 183'	Single Channel	Pad	10 Pads
32017197	BDK 007966	8.438" x 183' 200sh/pd	Single Channel	Pad	10 Pads
2018511C	BDK 007958	50mm x 150'	Single Channel	Roll	10 Rolls
31334595	BDK 007983	8.465" x 183' 200sh/pd	Single Channel	Pad	10 Pads
32916653	BDK 716-0237-00	8.465" x 183' 200sh/pd	Single Channel	Pad	10 Pads
32017188	BDK 716-0234-00	8.5" x 183' 200sh/pd	Multi Channel	Pad	10 Pads
32017526	BDK 007988 (716-0239-00)	215mm x 137.5' 300sh/pd	Multi Channel	Pad	10 Pads
32017526	BDK 716-0239-00	215mm x 137.5' 300sh/pd	Multi Channel	Pad	10 Pads
2026689C	BDK 007985	210mm x 60m 200sh/pd	Multi Channel	Pad	10 Pads

Corometrics


Product #	Cross Ref.	Description	Channels	UOM	Per Case
30597226	CMS 4305	152mm x 47' with Labels	Fetal	Pad	40 Pads
30767589	CMS 4305 BAO	152mm x 47' 160pg	Fetal	Pad	40 Pads
2003454C	CMS 4483	140mm x 130'	Fetal	Pad	12 Pads
4001556C	CMS 4305DAO 25 yr	152mm x 47'	Fetal	Pad	40 Pads

Datascope

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30015682	0683-00-0422-02 (SPA 307438)	50mm x 87'	Single Channel	Roll	50 Rolls

Dinamap

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018145	DM 89100	2.25" x 80' No Print	n/a	Roll	10 Rolls



Edan

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4007378C	EDN SE-12/1200	210mm x 295mm	Multi Channel	Pad	1 Pad

Electrodyne

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017187	ELD K140-1000	50mm x 89' Black Grid	Single Channel	Roll	100 Rolls
10669183	ELD K140-1010	50mm x 89' Red Grid	Single Channel	Roll	50 Rolls
10784107	ELD K140-1020	50mm x 89' Green Grid	Single Channel	Roll	50 Rolls

Fukuda Denshi

Product #	Cross Ref.	Description	Channels	UOM	Per Case
31052726	FKD 050F30-H02	50mm x 98'	Single Channel	Pad	120 Pads
2016931C	FKD OP-119TE	63mm x 97'	Single Channel	Roll	20 Rolls
4008207C	FKD OP-219TE	63mm x 99' 302sh/pd	Multi Channel	Pad	12 Pads
2025248C	FKD OP-124TE	49.5mm x 67' 272sh/pd	Multi Channel	Pad	10 Pads
2001532C	FKD OP-69TE	209mm x 88'	Multi Channel	Roll	10 Rolls

GE

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000936C	GE 2030887-001	80mm x 82' 280sh/pd	Single Channel	Pad	10 Pads
32018479	GE 74205-HEL	49.5mm x 84' No Print	n/a	Roll	20 Rolls
4001144C	GE 2036970-001 MAC 800	110mm x28m 200sh/pad	Multi Channel	Pad	10 Pads
4000936C	GE 2030887-001 MAC 400/600	79.5mm x 25.2m	Multi Channel	Pad	10 Pads

Getinge

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018706	GTG 61301600005	2.25" x 107' (can order 8)	n/a	Roll	80 Rolls
31449534	MRC GTG 61301601272	2.36 x 70' No Print	n/a	Roll	10 Rolls

Hellige

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2022033C	HEB 22616605	210mm x 45m 150sh/pd	Multi Channel	Pad	10 Pads
4000468C	HEB 22616611 25 yr	210mm x 45m 25yr	Multi Channel	Pad	1 Pad

Hewlett Packard

Product #	Cross Ref.	Description	Channels	UOM	Per Case
10005156	HP 9270-0484	151.6mm x 50'	Fetal	Pad	40 Pads
10632900	HP 40457A (can order 10)	50mm x 100'	Single Channel	Roll	50 Rolls
2000188C	HP 9270-0485	151.6mm x 50'	Fetal	Pad	40 Pads
20060622	HP 9270-0980 (can order 1)	50mm x 100'	Single Channel	Roll	50 Rolls
2000036C	HP 40477A	50mm x 89' No Print	Single Channel	Roll	10 Rolls
30748696	HP M1910A	5.894" x 49'	Fetal	Pad	40 Pads
30768678	HP M1707A	8.49" x 183' 200sh/pd	Multi Channel	Pad	8 Pads
4006120C	HP M1709A	209.5mm x 197'	Multi Channel	Pad	8 Pads
31104162	HP M2481A	8.484" x 183'	Single Channel	Pad	8 Pads
32016831	HP M1911A	150mm x 49'	Fetal	Pad	40 Pads
32021183	HP M1911A 25yr	150mm x 49' 25yr	Fetal	Pad	40 Pads

Hewlett Packard (Cont.)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018398	HP M1914A	150mm x 49'	Fetal	Pad	60 Pads
4000700C	HP 3708A	210mm x 30m 100sh/pd	Multi Channel	Pad	1 Pad

Humphrey

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2010895C	HUM 24433	216mm x 30m	n/a	Roll	1 Roll

Instrumentation Labs

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018666	ILS GEM 5508	59.5mm x 90' No Print	n/a	Roll	5 Rolls

Marquette | GE

Product #	Cross Ref.	Description	Channels	UOM	Per Case
01138221	MQE 9402-018	4.232" x 150'	Single Channel	Roll	32 Rolls
2000060C	MQE 9402-020	8.441" x 275' 300sh/pd	Multi Channel	Pad	1 Pad
2000151C	MQE 9402-021	8.234" x 275' 300sh/pd	Multi Channel	Pad	1 Pad
4006477C	MQE 9402-023	107.5mm x 75'	Multi Channel	Roll	2 Rolls
4004771C	MQE 9402-024	8.441" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
31002226	MQE 9402-024 GRN	8.443" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
31002200	MQE 9402-029	8.443" x 275'	Single Channel	Pad	8 Pads
31019568	MQE 9402-30	8.443" x 275' 300sh/pd	Multi Channel	Pad	8 Pads
32017205	MQE 9402-034 ARCHIVAL	8.433" x 275'	Multi Channel	Pad	12 Pads
2022282C	MQE 9402-037 ARCHIVAL	8.433" x 275'	Multi Channel	Pad	8 Pads
30781747	MQE 9402-046 (can order 10)	50mm x 89' No Print	Single Channel	Roll	50 Rolls
32017196	MQE 9402-053	8.5" x 11' 500sh/pd	Multi Channel	Pad	5 Pads
32016997	MQE 9402-061	8.433" x 139' 152sh/pd	Multi Channel	Pad	16 Pads

Medtronic

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32020245	MEN 6092	4.33" x 100' 204ah/pd	Multi Channel	Pad	6 Pads

Mennen

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30707080	MEG800-040-160	2.475" x 200'	Single Channel	Roll	50 Rolls

Micro Medical

Product #	Cross Ref.	Description	Channels	UOM	Per Case
0103931C	MI PSA1600	112mm x 8m	n/a	Roll	5 Rolls

Mindray

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4001129C	MRY 0683-00-0422-02 TR Paper	50mm x 18.75m	Single Channel	Roll	5 Rolls

Miscellaneous

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32016732	MED 705225 (LifeSigns 5200) (Can order 10)	57mm x 28'	Single Channel	Roll	50 Rolls

Mortara

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30758273	MRN 9100-007	8.465" x 458' 500sh/pd	Multi Channel	Pad	1 Pad
32020410	MRN 9100-025-50	215mm x 92' 100sh/pd	Multi Channel	Pad	10 Pads
4000482C	MRN 9100-026-50	8.433" x 62.72m 246/sh/pd	Multi Channel	Pad	1 Pad
0104415C	MRN 9100-026-51	8.250" x 75m 250sh/pd	Multi Channel	Pad	5 Pads
4003635C	MRN 9100-028-50	108mm x 140mm 200sh/pd	Multi Channel	Pad	15 Pads

Nihon Kohden

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000064C	NKN RQS50-3	50mm x 30m	Single Channel	Roll	20 Rolls
2010078C	NKN RQS63-3	63mm x 30m	Single Channel	Roll	5 Rolls
2026452C	NKN FQW110-2-140	110mm x 21m	Multi Channel	Pad	1 Pad
2026712C	NKN FQW210-30-140	210mm x 30m	Single Channel	Pad	10 Pads

Philips

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32018225	PHM 989803137011	8.464" x 11" 250sh/pd	Multi Channel	Pad	12 Pads
4000650C	PHM 989803138171	75mm x 25m	Single Channel	Roll	5 Rolls
4000651C	PHM 989803102321	215mm x 84m 300sh/pd	Multi Channel	Pad	5 Pads

Physio-Control

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2023274C	PHY 805319-03	106.5mm x 75'	Single Channel	Roll	2 Rolls
30718970	PHY 804700 (can order 5)	50mm x 100'	Single Channel	Roll	50 Rolls

Protocol

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32017249	PRT 008-0040-98 (Can order 40)	60mm x 87' No Print	Single Channel	Roll	100 Rolls

Quinton

Product #	Cross Ref.	Description	Channels	UOM	Per Case
0105967C	QUT 036868-001	210mm x 300mm 200sh/pd	Multi Channel	Pad	10 Pads
30561354	QUT 16211-002	8.5" x 183' 200sh/pd	Multi Channel	Pad	20 Pads
30729225	QUT 18416-001	8.5" x 183' 200sh/pd	Multi Channel	Pad	20 Pads
01083203	QUT 15022-001	8.484" x 200' 200sh/pd	Multi Channel	Pad	10 Pads
2011421C	QUT 30146-001	210mm x 60m 200sh/pd	Multi Channel	Pad	1 Pad

Radiometer

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2025244C	RMR 984-070	111.5mm x 142'	Single Channel	Roll	5 Rolls

Schiller

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000319C	SCH 2.157005	145mm x 115'	Single Channel	Pad	25 Pads
30999067	SCH 2.157007	70mm x 59' 200sh/pd	Single Channel	Pad	36 Pads
2000320C	SCH 2.157009	210mm x 100.8m	Multi Channel	Pad	12 Pads
2022430C	SCH 2.157014	90mm x 36m	Single Channel	Pad	1 Pad
2023563C	SCH 2.157016	90mm x 23m	Single Channel	Pad	12 Pads

Schiller (Cont.)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2023699C	SCH 2.157017	8.243" x 56m	Multi Channel	Pad	1 Pad
32016808	SCH 2.157023	210mm x 69' 150sh/pd	Multi Channel	Pad	20 Pads
0101533C	SCH 2.157025	8.243" x 48m	Multi Channel	Pad	10 Pads
0103082C	SCH 2.157026	80mm x 21m	Multi Channel	Pad	1 Pad

Siemens

Product #	Cross Ref.	Description	Channels	UOM	Per Case
20371540	SIS 43-55-060	60mm x 150'	Single Channel	Roll	50 Rolls
32018703	SIS 4398	4.37" x 80'	n/a	Roll	5 Rolls
4009491C	SIS 5773 (Clinitek 50/100/500)	2.24" x 79' No Print	n/a	Roll	5 Rolls
32019289	SIS 101404	2.250" x 150' No Print	n/a	Roll	10 Rolls

Sonicaid (Huntleigh)

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2000089C	SON 8400-8003	142.5mm x 45m	Fetal	Pad	1 Pad
4003692C	SON 8400-8003 25yr	142.5mm x 45m 25yr	Fetal	Pad	1 Pad

Spacelabs

Product #	Cross Ref.	Description	Channels	UOM	Per Case
30015682	SPA 307438 (Can order 10)	50mm x 90'	Single Channel	Roll	50 Rolls
30540499	SPA 307568	120mm x 150' 305sh/pd	Multi Channel	Pad	20 Pads
30763711	SPA 006-0196-00	50mm x 90'	Single Channel	Pad	100 Pads
4004322C	SPA 006-0210-02	57mm x 10m	n/a	Roll	12 Rolls

St. Johns

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000461C	SJM 60003605-001	210mm x 28m	n/a	Pad	12 Pads

Steris

Product #	Cross Ref.	Description	Channels	UOM	Per Case
32016258	MRC AMS P129359-008	2.28" x 75' No Print (can order 5)	n/a	Roll	100 Rolls
32008832	AMS P764330-673	American Sterilizer	n/a	Roll	24 Rolls

Welch Allyn

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4000944C	WA 94016-0000	8.243" x 56m 200sh/pd	Single Channel	Pad	5 Pads
4000505C	WA 94018-0000	8.243" x 56m 200sh/pd	Single Channel	Pad	1 Pads

Vitalograph

Product #	Cross Ref.	Description	Channels	UOM	Per Case
4007874C	VTO 44058	111mm x 35m	n/a	Roll	5 Rolls
30868940	VTO 20200	7.8125 x 8.594"	Single Channel	Pad	15 Pads

Zoll

Product #	Cross Ref.	Description	Channels	UOM	Per Case
2026560C	ZL 8000-0300	90mm x 18m 200sh/pd	Single Channel	Pad	1 Pad
4009107C	ZL 8000-000901 (X-Series)	79mm x 27.5m	Multi Channel	Roll	20 Rolls



Medical Chart Paper

Other Products We Offer

- **Sterilizer Charts:** Getinge & Steris
- **Lab Pharmacy:** Steris, Panasonic, Fujitsu, Linear Instruments, Pyxis, Seiko
- **Facilities Maintenance Charts:** Honeywell, Weksler, Bailey, United Electric, Instron
- **Variety of Pens and Markers**

Facilities Charts

Standardization Made Simple™

Helmer Blood Bank Refrigeration Charts

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0105872C	HER 220419 7 Day	4.5"	52	3 Boxes
0098882C	HER 220366 7 Day	4.5"	52	3 Boxes
0101104C	HER 220064 7 Day	4.5"	52	3 Boxes
0091359C	HER 220273 7 Day	4.5"	52	1 Box

Jewitt

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0099041C	JT RDR021 - Circle Chart	7.438"	52	3 Boxes
0097630C	JT RDR020 - Circle Chart	7.438"	52	3 Boxes
0098007C	JT RDR019 - Circle Chart	7.438"	52	3 Boxes
0097831C	JT RDR068 - Circle Chart	7.438"	52	3 Boxes
0097830C	JT RDR069 - Circle Chart	7.438"	52	3 Boxes


United Electric

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0084845C	UD 6282-169 - Circle Chart	6"	50	1 Box
0087720C	UD 6282-224 - Circle Chart	6"	50	3 Boxes
0096304C	UD 6282-176 - Circle Chart	6"	50	3 Boxes
0087703C	UD 6282-240 - Circle Chart	6"	50	3 Boxes
0087561C	UD 6282-161 - Circle Chart	6"	50	3 Boxes
0096344C	UD 6282-301 - Circle Chart	6"	50	3 Boxes


Jules Richard

Product #	Cross Ref.	Description	Charts / Box	Minimum Order
0075593C	RIC D32170 - Circle Chart	4.937"	100	1 Box
0084442C	RIC D32870 - Circle Chart	4.937"	100	1 Box
0100906C	RIC D33070 - Circle Chart	4.937"	100	1 Box
0078013C	RIC D31740 - Circle Chart	4.937"	100	1 Box
0089870C	RIC D32840 - Circle Chart	4.937"	100	3 Boxes
0101142C	RIC D32340 - Circle Chart	4.937"	100	3 Boxes

0084845C



0099041C



Ultrasound Gels



ULTRASOUND GELS

Vermed carries an extensive line of ultrasound gels and accessories. From the most widely used gel for diagnostic ultrasound - Aquasonic 100, to Gel Warmers and Disinfectants. We carry all the gels and accessories required throughout your facility.

Ultrasound Gels

Standardization Made Simple™

Vermed

Specially formulated to ensure excellent signal transmission

Product #	Cross Ref.	Description
4009235C (Pictured Right)	G250B	250ml Bottle - Blue ultrasound gel



Aquasonic 100

The most widely used gel for diagnostic and therapeutic medical ultrasound procedures

Product #	Cross Ref.	Description
4003889C (Pictured Right)	01-50	5L Sonicpac w/ Refillable Dispenser
4003891C	01-08	0.25L Dispenser
4003892C	01-02	60g Doppler sized tube
4007229C	01-34	1L w/ Dispenser and Cap



Aquasonic Clear

Contains no fragrance or colour. For diagnostic and therapeutic medical ultrasound procedures

Product #	Cross Ref.	Description
4000516C (Pictured Right)	03-50	5L Sonicpac w/ Refillable Dispenser
4008448C	03-08	0.25L Dispenser
4009439C	03-02	60g Doppler sized tube



4000516C

4003889C

4009235C

SCAN Ultrasound Gels

Less viscous gel


Product #	Cross Ref.	Description
4009440C	11-08	237ml (8 fl oz.) Dispenser
4009329C (Pictured Right)	11-28	3.8L (1 Gallon) w/ Refillable Dispenser



Ultrasound Probe Covers

Fits most endocavity ultrasound probes

Product #	Cross Ref.	Description
4006729C	38-01	Latex Free Probe Cover
32021026	38-03	3D/4D Probe Cover



Thermasonic Gel Warmers

Single and 3 Bottle Warmers

Product #	Cross Ref.	Description
32021021	82-01-20 CE	Single Bottle Warmer - 230V (Euro Plug)
32021022	82-03-20 CE	Multi Bottle Warmer - 230V (Euro Plug)
32024577 (Pictured Right)	82-03-20 CE	Multi Bottle Warmer LCD - 230V (Euro Plug)



Conductive Media

Electromedical Gels, Pastes and Adhesives

Product #	Cross Ref.	Description
4009203C (Pictured Right)	12-02	Spectra 360 Electrode Gel - 60g (2oz tube)
2000764C (Pictured Right)	12-08	Spectra 360 Electrode Gel - 250g (8.5oz tube)
32021031	22-06	Tensive Conductive Adhesive - 50g tube
32020230	67-05	Redux Electrolytic Paste - 5oz tube
2026134C	14-25	Signagel Electrode Gel - 250g (8.5oz tube)



Foetal Belts

Standardization Made Simple™

CTG / Abdominal Transducer Belts

Disposable and reusable foetal belts for secure CTG monitoring.

Product #	Cross Ref.	Description	Per Box
4009043C	DISPBELTBUT02	Disposable Abdominal Transducer Belt, 2 Pack - 1 x Pink, 1 x Blue, with 2 Buttons	100 Belts and 100 Buttons / 50 Packaged Pairs
4001531C	ROLLCTGBELT02	Disposable Abdominal Transducer Belt on a roll 52mm x 15m Black	-
4009543C	CTGBELT02	Reusable Abdominal Transducer Belt, 2 Pack Black	200 Belts / Box 100 Pairs
4004877C	DISPCTGBUT100	Disposable Abdominal Transducer Buttons	100 Buttons / Box

FOETAL BELTS

Designed for ease of use for clinicians and comfort for patients, these Transducer Belts have easy lock/unlock button holes, placed at intervals to obtain the most comfortable and safely secured CTG monitor position.

Latex-free



4001531C

4009543C

4009043C



Ultrasound Films

ULTRASOUND FILMS

NMT is an authorized distributor for Sony and Mitsubishi diagnostic imaging media, providing exceptional image detail, excellent resolution, and longevity. We also offer Vermed Ultra 110HD and Ultra 110HG as competitive alternatives compatible with Sony printers.

Ultrasound Films

Standardization Made Simple™

SONY

Sony ultrasound paper provides exceptional image detail, excellent resolution and longevity of the printed image. All of these outstanding qualities enable Sony to offer the most trusted products to explore what is happening inside the body with great precision.



MOST POPULAR

UP110 Series BLACK & WHITE



The best selling **UPP-110 series** offers the following advanced features: high water resistance, head-matching performance, gray-scale reproduction, anti-electrostatic layer, minimal curling, high humidity and heat resistance, superior print quality and excellent tearing properties.

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2021977C	UPP110HA	Superior Density Printing Paper FOR USE IN: UP-880, UP-890MD	110mm x 18m	240 Prints	10 Rolls
2019086C	UPP110HD	High Density Printing Paper FOR USE IN: UP-890MD, UP-895MD, UP-D895, UP-D898MD, UP-X898MD	110mm x 20m	270 Prints	10 Rolls
0098226C	UPP110HG	High Gloss Media FOR USE IN: UP-895MD, UP-D895, UP-D898MD, UP-X898MD	110mm x 18m	240 Prints	10 Rolls
2019085C	UPP110S/5	Standard Media FOR USE IN: UP-890MD, UP-895MD, UP-D895, UP-897MD, UP-D897, UP-D898MD, UP-X898MD	110mm x 20m	270 Prints	10 Rolls

SONY - Black and White Ultrasound Films

Product #	Cross Ref.	Description
2026054C	UPT210BL	Blue Thermal Transparency Film FOR USE IN: UP-980, and UP-991AD
0104937C	UPP725	8 x 10 Print Paper FOR USE IN: UP-D72XR
0106885C	UPT735BL	8 x 10 Blue Transparency Film FOR USE IN: UP-D71XR
31422960	UPP210SE	Standard Black & White Media FOR USE IN: UP-910, UP-930, UP-960, UP-980, UP-971AD, and UP-991AD
2026053C	UPP210HD	High Density Black & White Media FOR USE IN: UP-910, UP-930, UP-960, UP-980, UP-971AD, and UP-991AD

SONY - Colour Ultrasound Films

Product #	Cross Ref.	Description
32017149	UPC510C	A5 Colour Print Pack with Disposable Ribbon Holder FOR USE IN: UP-50, UP-D50, and UP-51MD
0104325C	UPC21L	A6 Colour Print Pack FOR USE IN: UP-21MD, UPD23MD, UP20, and UP25
0106244C	UPC21S	A6 Colour Print Pack FOR USE IN: UP-21MD, and UP-D21MD
32017030	UPC510	A5 Colour Print Pack FOR USE IN: UP-50, UP-D50, and P-51MD
32017150	UPC55	A5 Colour Print Pack FOR USE IN: UP-D55, and UP-D55MD

ULTRASOUND FILMS

MITSUBISHI

Mitsubishi ultrasound paper is a high quality product providing good resolution for ultrasound diagnostic printing. This paper provides great image detail when used with Mitsubishi printers such as: CP-700DSU, CP-700DU, CP-700UM, CP-770DW, P-68U, P-90W, P-91DW, P93DW, P-93W and P95DW.



7G30751336

Mitsubishi - Monochrome Ultrasound Films

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2018832C	K-61B	Super Grade Thermal Paper FOR USE IN: P-40U, P-50U, P-63UM, P-65U, P-67UA, P-68U, P-90U, P-91DW, P-91W, P93W, and P-95DW Mitsubishi Printers.	4.13" x 5.83"	100 Prints	20 Rolls
30957073	KP-65HM	High Density Thermal Sensitive Paper - Blue FOR USE IN: CP-700DSU, CP-700U, CP-700UM, CP-770DW, P-68U, P-90W, P-91DW, P-91W, P-93DW, P-93W, and P-95DW Mitsubishi printers.	110mm x 20m	215 Prints	20 Rolls
2020731C	K-75HM	High Density Thermal Paper - Blue FOR USE IN: P-71U, P-78U, and P-500W Mitsubishi printers	220mm x 25m	115 Prints	20 Rolls
0106115C	K-91HG	High Gloss Thermal Paper FOR USE IN: P-91DW, P93DW, P-93W and P-95DW Mitsubishi Printers.	110mm x 20m		20 Rolls

Mitsubishi - Colour Ultrasound Films

Product #	Cross Ref.	Description	Size	Per Pack	Per Case
2023136C	CK-700S	-	-	-	-
32017137	CK-700	Colour Roll Paper. 10mm x 22m FOR USE IN: CP-700DSU, CP-700DU, CP-700UM, CP-710U, and CP-770DW Colour printers	4" x 6"	-	10 Rolls
4008756C	CK-700L4P	Colour 4 Panel Paper Roll and Ink Ribbon FOR USE IN: CP-700DSU, CP-700DU, CP-700UM, CP-710U, CP-750U, and CP-770DW Colour printers	5.16 x 3.75"	75 Sheets	10 Rolls
0105985C	CK-3020L4P	-	-	-	-
2023137C	CK-700L	3 Panel Colour Paper Roll and Ink Sheet FOR USE IN: CP-800UM Colour printers	4 1/2 x 5"	160 Prints	10 Packs
4006030C	CK-800L	3 Panel Colour Paper Roll and Ink Sheet FOR USE IN: CP-800UM Mitsubishi printers	5" x 7"	100 Prints	10 Packs
0106262C	CK-900L	Standard L Size Ink Sheet FOR USE IN: CP-900DW, CP-900UM, and CP-910U Colour printers.	4.3" x 6.2"	130 Prints	10 Packs
0106263C	CK-900S4P	-	-	-	-
0104170C	CK-900S	Standard S Size Ink Sheet FOR USE IN: CP-900DW, CP900UM, and CP910U Colour printers.	4.1" x 4.3"	200 Sheets	10 Packs
4007230C	CK-3812	-	-	-	-
2019339C	CK-50S	-	-	-	-

ULTRASOUND FILMS

VERMED ULTRA

Available in both high density and high gloss formats, Vermed Ultra 110HD and Ultra 110HG are competitively priced thermal ultrasound films compatible with Sony medical grade printers.

Ultra 110HD - High Density Ultrasound Film

Product #	Rolls / Box	Rolls / Case
32019065	5 Rolls	50 Rolls

Ultra 110HG - High Gloss Ultrasound Film

Product #	Rolls / Box	Rolls / Case
32019066	5 Rolls	50 Rolls





European Headquarters

Torbay Business Park
Woodview Rd, Paignton
TQ4 7HP
United Kingdom
+44 1803 860100
UK.NisshaMedical.com

NORTH AMERICA • BUFFALO, NY | SANDPOINT, ID | COLORADO SPRINGS, CO | WOLCOTT, CT | GANANOQUE, ONTARIO

CARRIBEAN • DOMINICAN REPUBLIC

EUROPE • TORBAY, UNITED KINGDOM | DÜREN, GERMANY | RUM, AUSTRIA | GENT, BELGIUM | CERGY, FRANCE