



Defibrillator Compatibility

Match Defibrillator Brands to NMT-HS Part Numbers

BRAND OF DEFIBRILLATOR	MODEL OF DEFIBRILLATOR	NMT-HS COMPATIBLE DEFIBRILLATION ELECTRODE PART NUMBERS
Philips Medical	<ul style="list-style-type: none"> HeartStart Intrepid HeartStart MRx (Plug Style Connector) 	<ul style="list-style-type: none"> HeartStart XL T100AC-PHILIPS T100LOAC-PHILIPS C100AC-PHILIPS C100LOAC-PHILIPS PEDIATRIC-PHILIPS*
	<ul style="list-style-type: none"> HeartStart FR2+ DFM100 	<ul style="list-style-type: none"> HeartStart Forerunner T100AC-PHILIPS T100LOAC-PHILIPS C100AC-PHILIPS C100LOAC-PHILIPS
	<ul style="list-style-type: none"> HeartStart FRx 	<ul style="list-style-type: none"> HeartStart FR3 NOT COMPATIBLE
Physio-Control	<ul style="list-style-type: none"> LP-CR plus series LP-15 series 	<ul style="list-style-type: none"> LP-20 series LP-12 series T100AC-PHYSIO T100LOAC-PHYSIO C100AC-PHYSIO C100LOAC-PHYSIO PEDIATRIC-PHYSIO*
	<ul style="list-style-type: none"> LP-500 series 	<ul style="list-style-type: none"> LP-1000 series T100AC-PHYSIO T100LOAC-PHYSIO C100AC-PHYSIO C100LOAC-PHYSIO
Zoll Medical	<ul style="list-style-type: none"> M-Series E-Series R-Series 	<ul style="list-style-type: none"> Powerheart G3 Propaq MD X-Series T100AC-ZOLL T100LOAC-ZOLL C100AC-ZOLL C100LOAC-ZOLL PEDIATRIC-ZOLL*
	<ul style="list-style-type: none"> Zoll AED PLUS Zoll AED Pro 	<ul style="list-style-type: none"> Zoll AED 3 Zoll AED 3 BLS NOT COMPATIBLE
Cardiac Science	<ul style="list-style-type: none"> PowerHeart AED G3 Plus PowerHeart AED G3 Pro 	<ul style="list-style-type: none"> T100-CS
	<ul style="list-style-type: none"> Powerheart AED G5 	<ul style="list-style-type: none"> NOT COMPATIBLE
Schiller	<ul style="list-style-type: none"> FRED easy DG4000 DG5000 DG6002 	<ul style="list-style-type: none"> SKITY OBLS300 OBLS300A T100AC-SCHILLER T100LOAC-SCHILLER
	<ul style="list-style-type: none"> Fred Easy Fully and Semi Automatic AED 	<ul style="list-style-type: none"> Fred Easyport Pocket Defibrillator NOT COMPATIBLE
Nihon Kohden	<ul style="list-style-type: none"> TEC-5500 series TEC-7700 series TEC-5600 series 	<ul style="list-style-type: none"> TEC-8300 series TEC-7600 series T100AC-NIHONKOHDEN T100LOAC-NIHONKOHDEN
	<ul style="list-style-type: none"> CardioLife AED CardioLife AED-3100 	<ul style="list-style-type: none"> CardioLife AED Trainer CardioLife AED Linkage NOT COMPATIBLE

*HeartSync pediatric defibrillation electrodes are NOT compatible with AEDs. AEDs cannot be manually adjusted to reduce energy delivery. Please refer to your AED's manufacturer operation manual.

Please contact Nissha Medical Technologies-Healthcare Solutions directly for any questions regarding electrode compatibility.

© 2022 Nissha Medical Technologies. All Rights Reserved.

Disclaimer: The information presented is believed to be reliable concerning the subject matter covered. The final determination of the accuracy or completeness of any information, and the suitability of the product for the use contemplated, is the responsibility of the user.