



PAEDIATRIC & NEONATAL SINGLE-PATIENT NIBP CUFFS & HOSES

Nissha Medical Technologies provides a range of Non-Invasive Blood Pressure (NIBP) Cuffs and extensions especially made for Neonatal use. Our cuffs & hoses adapt to all brands and models of blood pressure monitors on the market. With more than 30 years of experience, NMT is the European expert in manufacturing and distribution of accessories linked to patient monitoring.



COMPLETE RANGE OF SIZES

Suits all types of patients from extra small premature to pediatric standard



FULL OEM COMPATIBILITY

Used with all brands on the market



QUALITY MATERIAL

Provides patient comfort



EXTENSIONS AVAILABLE

For all connectors on the market, allows the product to be ready for use in any circumstance



COLOUR CODED

Provides ease of identification and use



PAEDIATRIC/NEONATAL 1 TUBE CUFFS



Part Number	Patient	Cuff Color	Size	Per Box
98230305005	Paediatric	Orange	8-15cm	24
98230304005	Paediatric	Red	7-13cm	24
98230303005	Paediatric	Green	6-11cm	24
98230302005	Neonatal	Blue	4-8cm	24
98230301005	Neonatal	Pink	3-6cm	24

PAEDIATRIC/NEONATAL 2 TUBES CUFFS



Part Number	Patient	Cuff Color	Size	Per Box
98232305005	Paediatric	Orange	8-15cm	24
98232304005	Paediatric	Red	7-13cm	24
98232303005	Paediatric	Green	6-11cm	24
98232302005	Neonatal	Blue	4-8cm	24
98232301005	Neonatal	Pink	3-6cm	24



Provides consistent performance, with stable airtightness



Compliant with the latest regulations



A single-use product eliminating costly disinfection procedures and avoiding cross-contamination concerns



Transparent silk nylon material maximizing safety and providing patient comfort

NIBP HOSES

Nissha Medical Technologies offers a wide range of NIBP hoses that cover most of monitor brands on the market. The hoses dedicated to neonatal use have a pink color, differentiating them from the blue adult hoses.



Philips
Mindray



Nihon
Kohden



GE
Datex-Ohmeda



GE
Marquette



GE
Critikon