

ECG ELECTRODES



SWAROMED ELECTRODES


NISSHA
MEDICAL TECHNOLOGIES












FEATURES

Swaromed ECG Electrodes are pre-gelled (single-use) and have a Silver/Silver Chloride sensor. All adhesive foams are equipped with a hypoallergenic medical-grade adhesive that is non-irritating to the skin.



















- Hypoallergenic adhesive
- Silver/Silver Chloride sensor
- Quick recovery gel
- Strict quality control during production
- Defibrillation resistant










60mm - Round - Foam Electrodes

| Product # | Reference # | Per Pouch | Per Box | Per Case | Snap | Adhesive | Gel Type | | |
|--|-------------|-----------|---------|----------|-----------------|-----------|-----------|---|--|
|  6000003C | 1001 | 25 | 250 | 3000 | Stainless Steel | Regular | Wet Gel |  |  72 Hours |
| 6001122C | 4001 | 25 | 250 | 3000 | Stainless Steel | Regular | Solid Gel |  |  72 Hours |
| 6000004C | 1002 | 25 | 250 | 3000 | Stainless Steel | High Tack | Wet Gel |  |  48 Hours |
| 6001446C | 4002 | 25 | 250 | 3000 | Stainless Steel | High Tack | Solid Gel |  |  48 Hours |


















50mm - Round - Foam Electrodes

| Product # | Reference # | Per Pouch | Per Box | Per Case | Snap | Adhesive | Gel Type | | |
|--|-------------|-----------|---------|----------|-----------------|-----------|-----------|---|--|
|  6000018C | 1009 | 50 | 400 | 4800 | Stainless Steel | Regular | Wet Gel |  |  72 Hours |
| 6000582C | 4009 | 50 | 400 | 4800 | Stainless Steel | Regular | Solid Gel |  |  72 Hours |
| 6000436C | 1010 | 50 | 400 | 4800 | PP Carbon Blend | Regular | Wet Gel |   |  72 Hours |
| 6000981C | 4010 | 50 | 400 | 4800 | PP Carbon Blend | Regular | Solid Gel |   |  72 Hours |
| 6000665C | 1012 | 50 | 400 | 4800 | PP Carbon Blend | High Tack | Wet Gel |   |  48 Hours |
| 6000341C | 1019 | 50 | 400 | 4800 | Stainless Steel | High Tack | Wet Gel |  |  48 Hours |
| 6000616C | 4019 | 50 | 400 | 4800 | Stainless Steel | High Tack | Solid Gel |  |  48 Hours |

















55mm - Round - Foam Electrodes

| Product # | Reference # | Per Pouch | Per Box | Per Case | Snap | Adhesive | Gel Type | | |
|--|-------------|-----------|---------|----------|-----------------|-----------|-----------|---|--|
|  6000059C | 1023 | 50 | 400 | 4800 | Stainless Steel | Regular | Wet Gel |  |  72 Hours |
| 6000397C | 1022 | 50 | 400 | 4800 | Stainless Steel | High Tack | Wet Gel |  |  48 Hours |
| 6000652C | 4022 | 50 | 400 | 4800 | Stainless Steel | High Tack | Solid Gel |  |  48 Hours |

36 x 50mm - Oval - Foam Electrodes

| Product # | Reference # | Per Pouch | Per Box | Per Case | Snap | Adhesive | Gel Type | | |
|---|-------------|-----------|---------|----------|----------------------------------|-----------|-----------|---|--|
|  6000045C | 1036 | 50 | 500 | 6000 | Stainless Steel | High Tack | Wet Gel |  |  48 Hours |
| 6000578C | 4036 | 50 | 500 | 6000 | Stainless Steel | High Tack | Solid Gel |  |  48 Hours |
| 6000398C | 1040 | 50 | 500 | 6000 | PP Carbon Blend | High Tack | Wet Gel |   |  48 Hours |
| 6001237C | 4040 | 50 | 500 | 6000 | PP Carbon Blend | High Tack | Solid Gel |   |  48 Hours |
| 6000609C | 1084 | 30 | 300 | 3600 | Stainless Steel + Banana Plug | High Tack | Wet Gel |  |   48 Hours |
| 6000438C | 1084FT | 30 | 300 | 3600 | Stainless Steel + Banana Plug | High Tack | Wet Gel |  |   48 Hours |

32mm - Square & Teardrop - Foam Electrodes

| Product # | Reference # | Per Pouch | Per Box | Per Case | Snap | Adhesive | Gel Type | | |
|---|-------------|-----------|---------|----------|----------------------------------|-----------|-----------|---|--|
|  6000612C | 1057 | 30 | 300 | 3600 | Stainless Steel | Regular | Wet Gel |  |  72 Hours |
| 6000044C | 1060 | 50 | 500 | 6000 | PP Carbon Blend | High Tack | Wet Gel |   |  48 Hours |
| 6000588C | 4060 | 50 | 500 | 6000 | PP Carbon Blend | High Tack | Solid Gel |   |  48 Hours |
| 6000065C | 1066 | 50 | 500 | 6000 | Stainless Steel | High Tack | Wet Gel |  |  48 Hours |
| 6000581C | 4066 | 50 | 500 | 6000 | Stainless Steel | High Tack | Solid Gel |  |  48 Hours |
| 6000570C | 1083 | 50 | 400 | 4800 | Stainless Steel + Banana Plug | High Tack | Wet Gel |  |   48 Hours |

Banana Plug



Latex-free



Wet Gel



Solid Gel



MR Conditional



NISSHA

MEDICAL TECHNOLOGIES

United Kingdom

Torbay Business Park
Woodview Road
Paignton, Devon
TQ4 7HP
Phone: +44 1803 860100
Fax: +44 1803 863838
HS.NisshaMedical.com/uk/

Corporate Office

400 Exchange Street
Buffalo, NY 14204
Phone: 1.800.669.6905
Fax: 1.888.658.4941
HS.NisshaMedical.com