

ECG ELECTRODES



VERSATRODE ELECTRODES

NISSHA
MEDICAL TECHNOLOGIES



STAINLESS SNAP
Stops corrosion, eliminating artifact in the tracing. Silver / Silver Chloride sensor.

SUBSTRATE
1mm foam backing provides strength without stretching.

WET GEL
Chloride content optimised for reliable traces.



Nissha Medical Technologies' VersaTrode™ electrodes are the ideal choice for all adult cardiac monitoring and diagnostic applications. They are designed to work in all procedures: Stress Test, Long-Term Monitoring, Short Term Monitoring or Resting ECG.

60mm - Round - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type		
6000937N	FW60	25	250	3000	Stainless Steel	Regular	Wet Gel		72 Hours
6000939N	HTFW60	25	250	3000	Stainless Steel	High Tack	Wet Gel		48 Hours

50mm - Round - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type		
6000936N	FW50	50	400	4800	Stainless Steel	Regular	Wet Gel		72 Hours
6000935N	FW50R	50	400	4800	PP Carbon Blend	Regular	Wet Gel		72 Hours
6000941N	FS50	50	400	4800	Stainless Steel	Regular	Solid Gel		72 Hours
6000950N	FS50R	50	400	4800	PP Carbon Blend	Regular	Solid Gel		72 Hours
6001491N	HTFS50	50	400	4800	Stainless Steel	High Tack	Solid Gel		48 Hours
6001493N	HTFW50	50	400	4800	Stainless Steel	High Tack	Wet Gel		48 Hours
6001492N	HTFW50R	50	400	4800	PP Carbon Blend	High Tack	Wet Gel		48 Hours

36 x 50mm - Oval - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type		
6000982N	HTFW365B	30	600	3600	Stainless Steel	High Tack	Wet Gel		48 Hours
6000940N	HTFS365	50	500	6000	Stainless Steel	High Tack	Solid Gel		48 Hours

32mm - Square - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type		
6000930N	HTFW323	50	500	6000	Stainless Steel	High Tack	Wet Gel		48 Hours
6001023N	HTFS323	50	500	6000	Stainless Steel	High Tack	Solid Gel		48 Hours
6001126N	HTFS323R	50	500	6000	PP Carbon Blend	High Tack	Solid Gel		48 Hours
6001200N	HTFW323B	50	400	4800	Stainless Steel	High Tack	Wet Gel		48 Hours
6001494N	HTFW323R	50	400	4800	PP Carbon Blend	High Tack	Wet Gel		48 Hours

32mm - Teardrop - Foam Electrodes



Product #	Reference #	Per Pouch	Per Box	Per Case	Snap	Adhesive	Gel Type		
6000931N	FW32P	50	500	6000	Stainless Steel	Regular	Wet Gel		72 Hours
6001437N	FS32P	50	500	6000	Stainless Steel	Regular	Solid Gel		72 Hours
6001409N	FS32PR	50	500	600	PP Carbon Blend	Regular	Solid Gel		72 Hours

CARBON SNAP

MR Conditional and Radiolucent with a Silver / Silver Chloride sensor.

LABEL

Circular dot label holds snap in place on the electrode.

ADHESIVE

Medical grade acrylic adhesive ensures electrodes stay in place.

SOLID GEL

Aids in the adhesion of electrode and does not require additional cleanup.

FEATURES

- Adult foam; wet and solid gel electrodes
- Hypoallergenic adhesive
- High Tack option offers stronger adhesive
- For use throughout the hospital

Banana Plug



Latex-free



Wet Gel



Solid Gel



MR Conditional



NISSHA

MEDICAL TECHNOLOGIES

United Kingdom

Torbay Business Park
Woodview Road
Paignton, Devon, UK
TQ4 7HP
Phone: +44 1803 860100
[HS.NisshaMedical.com/uk](https://www.HS.NisshaMedical.com/uk)

Corporate Office

400 Exchange Street
Buffalo, NY 14204
Phone: 1.800.669.6905
Fax: 1.888.658.4941
[HS.NisshaMedical.com](https://www.HS.NisshaMedical.com)