



ZBrand™  
All Purpose Electrodes

**FEATURES**



**VARIETY OF SHAPES AND SIZES**  
Available in multiple shapes and sizes for better nesting of electrodes in close spaces.



**WET GEL OR SOLID GEL**  
Both gel options available to meet your clinical preference.



**LIFT TABS**  
Lift tab allows for easy removal from liner.

ZBrand™ electrodes were designed for a general patient monitoring solution with quality performance.


**PACKAGING**

**RESEALABLE BULK POUCH**

Ziplock pouch contains 30 or more electrodes

**CASE**

FDA / GS1 Compliant Labels



**HYGIENIC, DISPOSABLE SINGLE-USE PACKS**

Reduces risk of cross contamination of unused electrodes



# ZBrand™

All Purpose Electrodes


The Vermed ZBrand™ Electrode line was developed to meet the requirements of the entire hospital by offering an economic, general patient monitoring solution with quality performance.

## Z05MC - 40 x 36 mm - Teardrop Cloth Electrode




Product #	Ref #	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502683V	Z05MC-3	3	3	1200	●		72 Hours
6502684V	Z05MC-5	5	5	1000	●		72 Hours
6502685V	Z05MC-50	5	50	1000	●		72 Hours

## Z05MCRT - 50 x 48 mm - Round Cloth Electrode



Product #	Ref #	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502686V	Z05MCRT-50	5	50	1000	●	⚠	72 Hours

## Z05 - 40 x 36 mm - Teardrop Foam Electrode




Product #	Ref #	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502692V	Z05SG-3	3	3	1200	●		72 Hours
6502694V	Z05SG-5	5	5	1000	●		72 Hours
6502693V	Z05SG-30	3	30	1200	●		72 Hours
6502695V	Z05SG-50	5	50	1000	●		72 Hours
6502696V	Z05SGRT-30	3	30	1200	●	⚠	72 Hours
6502697V	Z05SGRT-50	5	50	1000	●	⚠	72 Hours
6502681V	Z05-5	5	5	1000	●		72 Hours
6502682V	Z05-50	5	50	1000	●		72 Hours

## Z35 - 40 x 35 mm - Rectangle Foam Electrode



Product #	Ref #	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502721V	Z35-3	3	3	1200	●		72 Hours
6502722V	Z35-5	5	5	1000	●		72 Hours
6502723V	Z35-50	5	50	1000	●		72 Hours
6502724V	Z35RT-3	3	3	1200	●	⚠	72 Hours
6502725V	Z35RT-5	5	5	1000	●	⚠	72 Hours
6502726V	Z35RT-50	5	50	1000	●	⚠	72 Hours

## Z05MT - 40 x 36 mm - Teardrop Tape Electrode



Product #	Ref #	Per Strip	Per Pouch	Per Case	Gel	MR Conditional	Duration
6502687V	Z05MT-3	3	3	1200	●		72 Hours
6502688V	Z05MT-5	5	5	1000	●		72 Hours
6502689V	Z05MT-50	5	50	1000	●		72 Hours
6502690V	Z05MTRT-5	5	5	1000	●	⚠	72 Hours
6502691V	Z05MTRT-50	5	50	1000	●	⚠	72 Hours