

The 4020/03 Pleated HEPA filter protects the patients airways effectively from exogenous microbial loads, reducing the risk of extrinsic colonization and infection.

## FEATURES



### LINE OF DEFENSE

Reduce the risk of extrinsic colonization and infection.



### HYDROPHOBIC FILTER MEMBRANE

Cannot pass potentially contaminated fluids under normal conditions.



### INHIBIT MICROORGANISMS

High bacterial and viral filtration efficiency rates to reduce risk of cross contamination





### UNIVERSAL SIZING

22mm x 22mm connector ports



## Bacterial and Viral Filtration to reduce Cross Contamination

### FILTER PLACEMENTS IN A VENTILATOR SYSTEM



Filter can be used by itself or in line with other filters.

### PACKAGING

#### POUCH

1 non-sterile filter per bag

#### CASE

50 Bags / case

### 4020/03 PLEATED HEPA FILTER

Part #	Filter Type	Dimensions	Efficiency	Connection
32031705	PLEATED HEPA FILTER	81.5mm x 68.5mm	BFE 99.9999% VFE 99.9999%	22M-22F

