



The AG7178 expiratory filter provides an additional line of defense against bacteria, viruses, and allergens that often plague CPAP/BiPAP and Ventilator patients. This disposable HEPA filter is designed to clean the outgoing air from Ventilators, CPAP, and BiPAP machines. It provides safe and effective filtration, helping patients breathe more naturally.



FEATURES



LINE OF DEFENSE

Protects against viruses and bacteria



EFFECTIVE FILTRATION

Simple and safe design



LARGELY COMPATIBLE

Works with most models



UNIVERSAL SIZING

22mm x 22mm connector ports



Making the World Safer, Healthier, and More Productive

AG7178 FINAL BACTERIAL / VIRAL FILTER FILTRATION PERFORMANCE

Part #	UOM	Filter Type	Filtration @0.3 micronB	BFE	VFE
32031690	Each	Electrostatic / HEPA	99.97%+	99.999%+	99.99%+

PACKAGING

INDIVIDUALLY PACKED

1 non-sterile filter per bag

CASE

100 bags / case

CONTACT

FOR MORE INFORMATION & AVAILABILITY

info@vermed.com

1.800.669.6905

Or contact your local Sales Representative.

