































Symbols Glossary Definition

Symbol	Standard Reference	Standard Title	Symbol Title	Description
	ISO 15223-1 Reference no. 5.4.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Do not re-use	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	ISO 15223-1 Reference no. 5.2.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Do not re-sterilize	Indicates a medical device that is not to be re-sterilized
	ISO 15223-1 Reference no. 5.2.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Sterilized using ethylene oxide	Indicates a medical device that has been sterilized using ethylene oxide.
	ISO 15223-1 Reference no. 5.1.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	ISO 15223-1 Reference no. 5.1.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	ISO 15223-1 Reference no. 5.1.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Serial number	Indicates the manufacturer's serial number so that a specific medical device can be identified.
	ISO 15223-1 Reference no. 5.1.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Date of manufacture	Indicates the date when the medical device was manufactured.
	ISO 15223-1 Reference no. 5.1.1	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Manufacturer	Indicates the medical device manufacturer, as defined in EU Directive 90/385/EEC, 93/42/EEC and 98/79/EC
	ISO 15223-1 Reference no. 5.1.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Use-by-date	Indicates the date after which the medical device is not to be used.

Symbol	Standard Reference	Standard Title	Symbol Title	Description
	ISO 15223-1 Reference no. 5.2.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Do not use if package is damaged	Indicates a medical device that should not be used if the package has been damaged or opened.
	ISO 15223-1 Reference no. 5.4.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Contains or presence of natural rubber latex	Indicates the presence of natural rubber or dry natural rubber latex as a material of construction within the medical device or the packaging of a medical device.
	ISO 15223-1 Reference no. 5.4.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
	ISO 15223-1 Reference no. 5.3.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.
	ISO 15223-1 Reference no. 5.3.6	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Upper limit of temperature	Indicates the upper limit of temperature to which the medical device can be safely exposed.
	ISO 15223-1 Reference no. 5.3.5	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Lower limit of temperature	Indicates the lower limit of temperature to which the medical device can be safely exposed.
	ISO 15223-1 Reference no. 5.1.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Authorized representative in the European Community	Indicates the authorized representative in the European Community
	ISO 15223-1 Reference no. 5.3.2	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Keep away from sunlight	Indicates a medical device that needs protection from light sources
	ISO 15223-1 Reference no. 5.3.4	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Keep dry	Indicates a medical device that needs protection from moisture
	ISO 15223-1 Reference no. 5.4.3	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Consult instructions for use	Indicates the need for the user to consult the instructions for use
	ISO 15223-1 Reference no. 5.2.7	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Non-sterile	Indicates a medical device that has not been subjected to a sterilization process.

Symbol	Standard Reference	Standard Title	Symbol Title	Description
	ISO 15223-1 Reference no. 5.3.8	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied – Part 1: General requirements.	Storage humidity range	Indicates the range of humidity to which the medical device can be safely exposed
	BS EN 15986:2011	Symbol for use in the labelling of medical devices — Requirements for labelling of medical devices containing phthalates.	Contains or presence of phthalate: bis (2-ethylhexyl) phthalate (DEHP)	Medical device is derived from or manufactured from products containing phthalate: bis (2-ethylhexyl) phthalate (DEHP).
	21 CFR 801.109	Labeling; Prescription devices.	Prescription Use Only	Caution: Federal law (USA) restricts this device to sale by or on the order of a licensed healthcare practitioner.
	765/2008/EC 768/2008/EC MDD 93/42/EEC Articles 4,11,12,17, Annex II)	The requirements for accreditation and market surveillance relating to the marketing of products; Medical Device Directive.	CE marking Conformité Européene Notified Body Reference no. ###	Signifies European technical conformity.
	n/a	n/a	PVC-Free	Medical device does not contain PVC
	n/a	n/a	MR Safe	Completely non-metallic, non-conductive, and not RF reactive.
	n/a	n/a	MR Conditional	The object or device is safe under certain tested conditions, and those conditions should be enumerated on the product, its packaging or in the enclosed literature
	ASTM F2503	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.	Magnetic Resonance (MR) Unsafe	An item which poses unacceptable risks to the patient, medical staff or other persons within the MR environment.
	n/a	n/a	Latex-Free	Indicates that natural rubber latex was not used in the manufacturing of the product, its container, or its packaging.
	n/a	n/a	Medical Device	Indicates that the product is a medical device.