



Protect eyes and face against liquids and dust. A full face shield offers maximum protection on all sides of the face, protecting eyes, nose and mouth from contact with liquids and dust. LATITUDE™ is not impact resistant nor intended for use against open flames or sparks.

## FEATURES



### LIGHTWEIGHT & OPTICALLY CLEAR

Provides clear, unobstructed view. Lightweight and clear. Effective for use against liquids and dust. LATITUDE™ PLUS: includes Anti-fog coating.



### MAXIMUM PROTECTION

Complete coverage including the sides of the face preventing all contact with eyes, mouth, and nose.





### COMFORTABLE AND FUNCTIONAL

LATITUDE™ comes with a comfortable brow foam. It also easily goes over glasses or goggles. LATITUDE™ is Latex, Fiberglass and Distortion Free.

## LATITUDE - PROTECTIVE FACE SHIELD

Product #	Description	Anti Fog	Per Bag	Per Case
32031704	<b>LATITUDE™ PLUS</b> Protective Face Shield 3M™ Anti-Fog Face Shield Film, 9960FS, 6.8 mil	☑	1	400
32031693	<b>LATITUDE™</b> Protective Face Shield Polyethylene Terephthalate (PET) Film	☒	1	400

## LATITUDE™ SPECS & PACKAGING



3-13-01-074 (5/2020)