



SPECIALTY FEATURES

AQUASONIC 100

The most popular gel recommended by leading manufacturers of medical ultrasound equipment worldwide.



AQUASONIC 100 SINGLE USE PACKS

Overwrapped, single use, sterilized foil pouches. Also available in non-sterile packs.

AQUASONIC CLEAR

Hypoallergenic gel with no fragrance or color.



SCAN

Less viscous gel that is applied as a gel and scanned as a fluid.

Parker ultrasound gel is bacteriostatic, non-sensitizing and non-irritating. Gels are completely aqueous and will not stain clothing or damage equipment. They are also acoustically correct, providing superior transmission of ultrasound waves. Ultrasound gel is not a spermicide and does not contain formaldehyde.



12 bottles / case



1 bottle / case
with refillable dispenser & cap



1 bottle / case
with refillable dispenser & cap



100 non-sterile packets / box
4 boxes / case

Parker Ultrasound Gels

Standardization Made Simple™

AQUASONIC 100 ULTRASOUND GEL

World Standard for Medical Ultrasound Transmission, AQUASONIC 100 is the most widely used gel for diagnostic and therapeutic medical ultrasound. Recommended for all procedures where a viscous gel is required.

Product #	Cross Ref	Description	# / Box
7G30592052	01-50	5L Sonipac with refillable dispenser	1
7G30592037	01-08	0.25L Dispenser	12
7G30900156	01-02	60g Doppler size tube	12
7G30981271	01-34	1L with dispenser and cap	6
7G30936127	01-01	20g Sterile packet	48
7G32021003	01-20	20g Non-sterile single use packs	100



7G30592052

7G32021003

AQUASONIC CLEAR ULTRASOUND GEL

CLEAR ultrasound gel. Based on the same “can’t be copied” formula of the World Standard, Aquasonic® 100 Ultrasound Transmission Gel, without fragrance or color. For diagnostic and therapeutic medical ultrasound procedures.

Product #	Cross Ref	Description	# / Box
7G30997020	03-50	5L Sonipac with refillable dispenser	1
7G30997012	03-08	0.25L Dispenser	12
7G32021004	03-02	60g Doppler size tube	12
7G32021005	03-34	1L with dispenser and cap	6
7G32021007	03-54	5L Econopac (4-5L dispensers)	1



7G30997020

SCAN ULTRASOUND GEL

Parker’s exclusive dual-viscosity formula allows SCAN to be applied as a gel and scan as a fluid. For medical ultrasound procedures when a less viscous gel is desired.

Product #	Cross Ref	Description	# / Box
7G30592078	11-08	237ml (8 fl oz) Dispenser	12
7G30592086	11-28	3.8L (1 Gallon) with refillable dispenser	1
7G30900198	11-28S	SCANPAC - 4 gallons with 2 dispensers and 1 pump	1



7G30592078