



SWAROPLATE

Non Split Neutral Electrode
1203H-swaro-V5-50-600
Reference: **1203H**

5 Electrodes / Pouch
600 Electrodes / Case

BIOCOMPATIBILITY	
ISO 10993	Passed
Latex Free	Yes

ENVIRONMENTAL	
RoHS Compliant	Yes
REACH Compliant	Yes

SHELF LIFE	
Product Shelf Life (in accordance with storage guidelines)	24 Months

ELECTRICAL PERFORMANCE (IEC 60601-2-2)	
Impedance by 600 kHz	4 Ohm

MATERIALS	
Substrate Material	PE Foam
Adhesive	Conductive Hydrogel
Release Liner	PETP Film
Connector	N/A
Cable	N/A

APPLICATION	
Grounding Pad for High Frequency Surgery for infants and pre-matures < 5kg	

PRODUCT INFORMATION	
Size	6.50 x 2.91 in / 165 x 74 mm
Substrate Thickness (adapter excluded)	0.039 in / 1 mm
Total Product Surface Area	10.23 in ² / 66 cm ²
Conductive Area	5.43 in ² / 35 cm ²
Integrated Lead Wire (length \ color)	N/A

PACKAGING	
Product Packaging Material	Foil - PET/AL/PE
Resealable Pouch	No
Product Packaging Size (L x W)	10.24 x 6.30 in 26 x 16 cm
Department Packaging - Box (Size / Qty)	12.20 x 6.69 x 2.95 in 31.0 x 17.0 x 7.5 cm Per Box: 50 e.a.
Transport Packaging - Carton (Size / Qty)	24.80 x 18.70 x 7.09 in 63.0 x 47.5x 18.0 cm Per Carton: 600 e.a.

Vermed Part Number: 6000788N