



SWAROPLATE

Split Neutral Electrode

1223H-swaro-V5-50-600

Reference: **1223H**

5 Electrodes / Pouch
600 Electrodes / Case

BIOCOMPATIBILITY

ISO 10993 Passed
Latex Free Yes

ENVIRONMENTAL

RoHS Compliant Yes
REACH Compliant Yes

SHELF LIFE

Product Shelf Life (in accordance with storage guidelines) 24 Months

ELECTRICAL PERFORMANCE (IEC 60601-2-2)

Impedance by 600 kHz 4 Ohm

APPLICATION

Grounding Pad for High Frequency Surgery for infants and pre-matures < 5kg

PRODUCT INFORMATION

Size 6.50 x 2.91 in / 165 x 74 mm
Substrate Thickness (adapter excluded) 0.039 in / 1 mm
Total Product Surface Area 10.23 in² / 66 cm²
Conductive Area 4.96 in² / 32 cm²
Integrated Lead Wire (length \ color) N/A

MATERIALS

Substrate Material PE Foam
Adhesive Conductive Hydrogel
Release Liner PETP Film
Connector N/A
Cable N/A

PACKAGING

Product Packaging Material

Foil - PET/AL/PE

Resealable Pouch

No

Product Packaging Size (L x W)

10.24 x 6.30 in
26 x 16 cm

Department Packaging - Box (Size / Qty)

12.20 x 6.69 x 2.95 in
31.0 x 17.0 x 7.5 cm

Per Box: 50 e.a.

Transport Packaging - Carton (Size / Qty)

24.80 x 18.70 x 7.09 in
63.0 x 47.5 x 18.0 cm

Per Carton: 600 e.a.

Vermed Part Number: 6000804N