


SWAROPLATE

Split Neutral Electrode

1252H-swaro-V1-8-96

Reference: **1252H**

1 Electrodes / Pouch
96 Electrodes / Case



BIOCOMPATIBILITY

ISO 10993 Passed
Latex Free Yes

ENVIRONMENTAL

RoHS Compliant Yes
REACH Compliant Yes

SHELF LIFE

Product Shelf Life (in accordance with storage guidelines) 24 Months

ELECTRICAL PERFORMANCE (IEC 60601-2-2)

Impedance by 600 kHz 2 Ohm

APPLICATION

Grounding Pad for High Frequency Surgery for children > 5kg and < 15kg

PRODUCT INFORMATION

Size 6.06 x 3.46 in / 154 x 88 mm
Substrate Thickness (adapter excluded) 0.039 in / 1 mm
Total Product Surface Area 19.22 in² / 124 cm²
Conductive Area 11.01 in² / 71 cm²
Integrated Lead Wire (length \ color) 3m \ blue

MATERIALS

Substrate Material PE Foam
Adhesive Conductive Hydrogel
Release Liner PETP Film
Connector Valleylab-REM
Cable PVC

PACKAGING

Product Packaging Material

Foil - PET/AL/PE

Resealable Pouch

No

Product Packaging Size (L x W)

11.81 x 6.30 in
30 x 16 cm

Department Packaging - Box (Size / Qty)

12.20 x 6.69 x 2.95 in
31.0 x 17.0 x 7.5 cm

Per Box: 8 e.a.

Transport Packaging - Carton (Size / Qty)

24.80 x 18.70 x 7.09 in
63.0 x 47.5 x 18.0 cm

Per Carton: 96 e.a.

Vermed Part Number: 6000833N