


TWIN-SAFE

Split Neutral Electrode
1285H-swaro-V1-8-96
Reference: **1285H**

1 Electrodes / Pouch
96 Electrodes / Case



BIOCOMPATIBILITY

ISO 10993 Passed
Latex Free Yes

ENVIRONMENTAL

RoHS Compliant Yes
REACH Compliant Yes

SHELF LIFE

Product Shelf Life (in accordance with storage guidelines) 24 Months

APPLICATION

Grounding Pad for High Frequency Surgery for adults and children

ELECTRICAL PERFORMANCE (IEC 60601-2-2)

Impedance by 600 kHz 1 Ohm

PRODUCT INFORMATION

Size 6.93 x 4.80 in / 176 x 122 mm
Substrate Thickness (adapter excluded) 0.039 in / 1 mm
Total Product Surface Area 26.20 in² / 169 cm²
Conductive Area 17.05 in² / 110 cm²
Integrated Lead Wire (length \ color) 3m \ blue

MATERIALS

Substrate Material PE Foam
Adhesive Conductive Hydrogel
Release Liner PETP Film
Connector Valleylab-REM
Cable PVC

PACKAGING

Product Packaging Material

Foil - PET/AL/PE

Resealable Pouch

No

Product Packaging Size (L x W)

11.81 x 6.30 in
30 x 16 cm

Department Packaging - Box (Size / Qty)

12.20 x 6.69 x 2.95 in
31.0 x 17.0 x 7.5 cm
Per Box: 8 e.a.

Transport Packaging - Carton (Size / Qty)

24.80 x 18.70 x 7.09 in
63.0 x 47.5 x 18.0 cm
Per Carton: 96 e.a.

Vermed Part Number: 6000845N