


TWIN-SAFE

Split Neutral Electrode

1286H-swaro-V1-8-96

Reference: **1286H**

1 Electrodes / Pouch
96 Electrodes / Case



BIOCOMPATIBILITY

ISO 10993	Passed
Latex Free	Yes

ENVIRONMENTAL

RoHS Compliant	Yes
REACH Compliant	Yes

SHELF LIFE

Product Shelf Life (in accordance with storage guidelines)	24 Months
---	-----------

APPLICATION

Grounding Pad for High Frequency Surgery for adults and children

ELECTRICAL PERFORMANCE (IEC 60601-2-2)

Impedance by 600 kHz	1 Ohm
-----------------------------	-------

PRODUCT INFORMATION

Size	6.93 x 4.80 in / 176 x 122 mm
Substrate Thickness (adapter excluded)	0.039 in / 1 mm
Total Product Surface Area	26.20 in ² / 169 cm ²
Conductive Area	17.05 in ² / 110 cm ²
Integrated Lead Wire (length \ color)	3m \ blue

MATERIALS

Substrate Material	PE Foam
Adhesive	Conductive Hydrogel
Release Liner	PETP Film
Connector	jack
Cable	PVC

PACKAGING

Product Packaging Material

Foil - PET/AL/PE

Resealable Pouch

No

Product Packaging Size (L x W)

11.81 x 6.30 in
30 x 16 cm

Department Packaging - Box (Size / Qty)

12.20 x 6.69 x 2.95 in
31.0 x 17.0 x 7.5 cm

Per Box: 8 e.a.

Transport Packaging - Carton (Size / Qty)

24.80 x 18.70 x 7.09 in
63.0 x 47.5 x 18.0 cm

Per Carton: 96 e.a.

Vermed Part Number: 6000846N