


# TWIN-SAFE

Split Neutral Electrode

1288H-swaro-V1-8-96

Reference: **1288H**

1 Electrodes / Pouch  
96 Electrodes / Case



## BIOCOMPATIBILITY

**ISO 10993** Passed  
**Latex Free** Yes

## ENVIRONMENTAL

**RoHS Compliant** Yes  
**REACH Compliant** Yes

## SHELF LIFE

**Product Shelf Life** (in accordance with storage guidelines) 24 Months

## ELECTRICAL PERFORMANCE (IEC 60601-2-2)

**Impedance by 600 kHz** 4 Ohm

## APPLICATION

Grounding Pad for High Frequency Surgery for infants and pre-matures < 5kg

## PRODUCT INFORMATION

**Size** 7.13 x 2.99 in / 181 x 76 mm  
**Substrate Thickness (adapter excluded)** 0.039 in / 1 mm  
**Total Product Surface Area** 11.94 in<sup>2</sup> / 77 cm<sup>2</sup>  
**Conductive Area** 5.74 in<sup>2</sup> / 37 cm<sup>2</sup>  
**Integrated Lead Wire (length \ color)** 3m \ blue

## MATERIALS

**Substrate Material** PE Foam  
**Adhesive** Conductive Hydrogel  
**Release Liner** PETP Film  
**Connector** Valleylab-REM  
**Cable** PVC

## PACKAGING

### Product Packaging Material

Foil - PET/AL/PE

### Resealable Pouch

No

### Product Packaging Size (L x W)

11.81 x 6.30 in  
30 x 16 cm

### Department Packaging - Box (Size / Qty)

12.20 x 6.69 x 2.95 in  
31.0 x 17.0 x 7.5 cm

**Per Box:** 8 e.a.

### Transport Packaging - Carton (Size / Qty)

24.80 x 18.70 x 7.09 in  
63.0 x 47.5 x 18.0 cm

**Per Carton:** 96 e.a.

**Vermed Part Number: 6000847N**